

Product datasheet for **MG214675**

1700080E11Rik (NM_028562) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700080E11Rik (NM_028562) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: 1700080E11Rik
Synonyms: A1429110; Tseg8
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG214675 representing NM_028562
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCACGCCCACTAAGATAGACTTGGAGGACGCCCTGGCCCTGGGCCCAGCTGGCGGCATCTTTGTC
 ACGTGATGCTGTATGCGCCCGATCCAGGGGTGCTCTTTGGTGGCCGTATCTTGCTGAACATGCTGTGCT
 GATGTCCATGCGCTTTGATGGGCGACTGGGCTCCCTGGTGGTTTTGTGGACGATCGTAGCCCCAGCCTG
 GAGGAAGGGCTGAATAAGGAAGTGTGAAGAACTGGGCGAGGGAGTGCCACTTTCAGCATCTTATCCA
 CCGACTACCGAAGCTCCCTCGCTGAGTCCAATCAAAGTGGTGGCCATTTCTATGTCAAATGCCTAAC
 ACTAGAGCAGCTCCAGGCGGTGGAGGCTGGGGCACCACAAGCCAAGGACCATGGGCTAGAGGTCTTGGGC
 CTGGTGCGGGTACCCCTGTACACCTTGC GGATGGTGTGGGAGGCCTGCCTGCTTTTCTGGAGAATTCT
 TCATTGGTGTGCGCCAGGGAACAGCTACTGGAAGCCCTCCAGGACTTGGGAATTCTGGCTCCCGAACTGC
 CTCAGACCTCAAGGAACGAACAGTTCAACGTAACAGACCTTAAGGCTGATTGGGACC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG214675 representing NM_028562
 Red=Cloning site Green=Tags(s)

MATPTKIDLEDALALGPSWRHLCHVMLYAPDPGVLFGRILLNYAVLMSMRFDGRLGFPGGFVDDRSPSL
 EEGLNKELLKLGEGVSTFSILSTDYRSSLAESKSKVVAHFYVKCLTLEQLQAVEAGAPQAKDHGLEVLG
 LVRVPLYTLRDGVGGLPAFLNSFIGVAREQLLEALQDLGILAPETASDLKERTVQRKQTLRLIGT

TRTRPLE - GFP Tag - V


Restriction Sites: Sgfl-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_028562

ORF Size: 618 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028562.3](#), [NP_082838.2](#)

RefSeq Size: 871 bp

RefSeq ORF: 621 bp

Locus ID: 73532

Cytogenetics: 9 F1