

Product datasheet for **MG214725**

1700037C18Rik (NM_028484) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: TurboGFP

Symbol: 1700037C18Rik

Synonyms: 2810022L20Rik

Mammalian Cell Selection: Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG214725 representing NM_028484
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGATGCAGTGAGCTGGGCCCCAGGACGCCCCAGCCACCCTGACACACCACCAACATCTATGAGGGAG
GATTGGGGGCCAGCAGCAGCAGGGGCCAGTGCTCAGGGAAGCAAGCCCAAGAAGCTCCGACTGCGACA
CCTCCGGAGCTGGCTCTCTACCTGCCAGGCCACATGCAGCCTGCTGGCCAATGTGGAAGCCACTGGCTA
GGCCGGCTCATGGCTGGGGCAGCCTGACACGGCCTGAAGGCTCACCCCTGGCCCTGGACCTGCCACAGG
GGACCCTGGGTCCAGGTAACAGCCATTGCTCAGCCCTTCTGGAAGCCCACTTGCCAGGGACAGCCTGGG
AAATACAGCTTCCAGTTCAGCATGGACCCAGCCAAGGGTGTCCCCTCCAGTCTGGTCCCCCGAGGGC
TTGGGACTCAGGCCAAGAGGTCCTGGAGAGCCTTGGAGGAGACCATGTGTCCCTTGTGCAAGAGAACC
GCTCTGGGCCTTGGAGAGGACA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG214725 representing NM_028484
Red=Cloning site Green=Tags(s)

MDAVSWAPGRPSHPDTPPNIEGGLGAQQQQPSAQQSKPKNFRLRHLRSLALYLPGHMQPAGQCQGS
HWRGRLMAGGSLTRPEGSPWPLDLPQGTLPGNSHCSALLEAHLPRDSLGNASSSSMDPAKGVPSQSGP
PEGLGLRPKRSWRALEETMCPLCKRTRSGALERT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_028484

ORF Size: 513 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_028484.2](#), [NP_082760.2](#)

RefSeq Size: 1099 bp
 RefSeq ORF: 516 bp
 Locus ID: 73261
 UniProt ID: [Q8BT88](#)
 Cytogenetics: 16 A1

Product images:



Circular map for MG214725