

Product datasheet for **MG214751**

Fam166b (NM_001162381) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam166b (NM_001162381) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Fam166b
Synonyms:	4833436C18Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG214751 representing NM_001162381 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGTGGTCAACACTTCCATCCCAGGGCTCAGCGGCGAGAATCCGCATTATATCCAGGGTACACCG
GACTGCCCCGCTACTTCGGTTCAGCATGGCCAGACCTATGGGCAGGTGACTGGTCAGCTACTCCGAGG
CCCTCCTGGTCTAGCCTGGCCTCCTGCCACCGCACACTTCTGCCTCCATCCAGTCTCCAAGATCTCCT
GTGATCTCCAAGGGGAGGCTGCCTCCCAGGCGCGGCATGAAAGGCTCAGCTCCAGCATAATCCCTGGT
ATACAGGTTTTATACCACGGGCACAGTTCAATTTTTGCCAAGAAGTCAACCAGGTTTGGGCTGAGCTAT
GAGTGAATTTACACGAAGGCACGGGAGCAAGAGAGCCATCAACTGCCAGACGGAGCCAAGGGAGAGAGA
GAAGTGAAGAAGACCAGCTGCGGGAGGCAAGAGGCCACCGCTGAAGCAGGAGTTGGCAGATGCTTCCC
CCTACTCCATGGATGATACAGATCCCCACAAGTTCTTTCATGTCAGGTTTACGGGCTATGTGCCCGTGC
CCGCTTCTCTTTGGCTCCAGCTTCCCTGTGCTTACCAACCAGGCATTGCAGGAATTTGGGCAGATGTGC
TCACGGGGCCGGGCGACAAGGACCCAAACCTCTCCCCACTTCCCAGGCCACGTTTCAGAATCTGG
GTCTCTACCTCACTACGGAGGCTACGTGCCAGGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG214751 representing NM_001162381
 Red=Cloning site Green=Tags(s)

MAVVNTSIPGLSGENPHYIPGYTGHCPLLRF SMGQTYGQVTGQLLRGPPGLAWPPAHR TLLPPIQSPRSP
 VISKGRLLP RRGHERLSSSIIPGYTGFIPRAQFIFAKNCNQVWAEAMSEFTRRHGEQESHQLPDGAKGER
 EVEEDQLREAE EPLKQELAHASPYSMDDTDPHKFFMSGFTGYVPRARFLFGSSFPVLTNQALQEF GQMC
 SRGRAHKDPKPLSPLPRPTFQNLG LLLPHYGGYVPG

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001162381

ORF Size: 735 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001162381.1 , NP_001155853.1
RefSeq Size:	1271 bp
RefSeq ORF:	738 bp
Locus ID:	329831
UniProt ID:	A2AIP0
Cytogenetics:	4 A5