

Product datasheet for **MG214760**

Spint5 (NM_001040055) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: Spint5
Synonyms: Gm14316; OTTMUSG00000016293; RP23-370H21.10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG214760 representing NM_001040055
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGGTCCAGGGATTGCGCTTGCCCTGTGGCTTCTAACGTTTGCCATGTCCCTGTCAATGCTATCTA
GTAGTTTGCTATATCCTGCAGGAGTCCGAAGTCAACTCTGCGAAAGCGGGCACCTTGAGGCGAAAAGAAT
AATATGCAACCAACCGGTAAGAAGGGCTTCTGCAGCTTCACTTTCTACAGATACTACTTCAACCCAGAA
TCGGCCCTATGTGAGTCCTTCATCTTCACTGGCTGCGGGGCAACAGAAACAACCTTTAAACTAAATACC
TCTGTGAAGTCCGTTGCATTACGTAGAGAAGGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG214760 representing NM_001040055
Red=Cloning site Green=Tags(s)

MAGPGIRLALWLLTFAMLSMLSSLLYPAGVRSQLCESGHLGAKRIICNQPVKKGFCSTFYRYFNPES
SALCESFIFTGCGGNRNNFKTKYLCEVRCIHVEKD

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: NM_001040055

ORF Size: 315 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001040055.1](#), [NP_001035144.1](#)


RefSeq Size: 443 bp

RefSeq ORF: 318 bp

Locus ID: 641368

Cytogenetics: 2 H3

Product images:



Circular map for MG214760