

## Product datasheet for **MG214769**

### 4931428F04Rik (NM\_001166394) Mouse Tagged ORF Clone

#### Product data:

<b>Product Type:</b>	Expression Plasmids
<b>Tag:</b>	TurboGFP
<b>Symbol:</b>	4931428F04Rik
<b>Synonyms:</b>	AI426165; Kiaa0895l
<b>Mammalian Cell Selection:</b>	Neomycin
<b>Vector:</b>	pCMV6-AC-GFP (PS100010)
<b>E. coli Selection:</b>	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:** >MG214769 representing NM\_001166394  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGTGCTGGACTCGGGACACAGGTGTATGAGCAGGCACCCCCAGGCCCTGCGGTTCCCATCTC  
AGCACCATAAGCTGAAGCCCTCAAACGAAATGGGCCACCCCTGTACCCCTGGCCTGAGTCCCTGGGGAT  
GCCCTTGGCTCTGGCAGTCCCTTCAGCACTGCAGCAGCAGACTATGTGGCAGACCTTCTCAAAGTGCAC  
TTGGAGCAGAGTAGCCACATGCGTCGACGAGAGCACCTACAGTGTAAATAGTACTGGCCGGCGGGGGC  
GTGGTAAAGCCCCGATTGGACGGGGATGCGATCCAGGAGGGACCCACGGCCTGCAGCCTCCTTACCTCA  
CATTGCTAAGATTGAAAAGGATGTAGGCAGTAGTAGCAGCAAGAGCCCTGTATGTTGGTGGCCCTGCGG  
CCAACATAACATGGACCAGGAACGGGAAAAGTTCTTCCAGTCTCATTATACCTACAATCCACAGTTCGAGT  
ATCAGGAACCGATGCCAATGAGTGTGCTGGAGAAGTACCAGGAGGCCTCTGCGCAATTCATGAATCAGGC  
AGTTGGCATCATTGAGGCTGTCTGGAGAAGTTTGGCACTTATGAAAACTTTGGGCTGCCACAGGGGGC  
CAGCTGCTGACCAAGTGCCAAATTTGGTCCACCGTGCGCAAATACATGCAGAAGGAAGGCTGTGTTGGGG  
AGATTGTGGTACAGCTGAGTGAAGTCTACTGTCCCAGGCTGTGATGATGGTGGAGAACAGCCGACCCAC  
ACTGGCCATCAACTTGACTGGAGCCCGCCAGTACTGGTTGGAGGGCATGCTGCGGCACGAGATAGGCACC  
CACTACCTGCGTGGTGTGAATAATTCGCGCAACCGTGGCACAGCACTGAAGGCAGATTGCAGTACGGCC  
TTCAGCCAGCCAACCCCTACTGAGGAGGGCCTGGCTAGCCTGCATAGTGTGCTGTTTCGAAAAGCAACCATT  
CCTGTGGCGCGCCGCCCTGCTCTACTACACCATCCACCAGGCGGCACATGTCCTTCCGCCAGCTCTTC  
CAGGATCTGGCGCAGTACGTGCAGGACGAGGGTGTACGCTGGGAGTACTGCGTCCGTGCCAAGCAGGAC  
AGACGGATACCTCCAGCCAGGCTGCTTACGCAAGGACCAGGTGTACCTGGATGGCATTTCGCGCATCCT  
ACGGCACCGGCACACCATTGATTTCCAGCTGCTGACCTCGCTGGGCAAGGTCTCCTATGAAGACGTGAA  
CAACTGCGGCCCTACGGGGTGTAGACAATACCAGAGTGCCCCACTTCATGAAGGACTTGAACGTTACC  
GGCAGCAGCTGGAACACATCATGACCACCAACCGGCTGGATGAAGCGGAGCTGGGCCGCTGCTGCCTGA  
C

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


**Protein Sequence:** >MG214769 representing NM\_001166394  
Red=Cloning site Green=Tags(s)

MVLDSGTQVYEQAPRPPAGSPSQHHKLPKPSNNGPPLYPWPESLGMPALAVPSALQQQTMWQTFSLH  
LEQSSHMRRSESTYSVNSTGRRGRGKAPIGRGCDPGGTLRPAASLPHIAKIRKDVGSSSSKSPCMLVALR  
PTNMDQEREKFFQSHYTYNPQFEYQEPMPMSVLEKYQEASAQFMNQAVGIEAVLEKFGTYENFEATGG  
QLLTKCQIWSVTRKYMKEGCVGEIVVQLSEDLLSQAVMMVENSRLPLAINTGARQYWLEGLRHEIGT  
HYLRGVNNSRQPWHSTEGRLQYGLQPANPTEGLASLHSLVFRKQPFLWRAALLYTIHQAAHMSFRQLF  
QDLAQYVQDEGRWEYCVRAKRQGTDTSPGCFSKDQVYLDGILRLRHRHTIDFQLLSLKGVSIEDVE  
QLRYPGVLDNTRVPHFMKDLERYRQLEHIMTTNRLDEAELGRLLPD

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001166394

**ORF Size:** 1401 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_001166394.1](#), [NP\\_001159866.1](#)

RefSeq Size: 3174 bp  
 RefSeq ORF: 1404 bp  
 Locus ID: 74356  
 UniProt ID: [Q810A5](#)  
 Cytogenetics: 8 D3

**Product images:**



Circular map for MG214769