

Product datasheet for **MG214800**

Zfp651 (NM_001166644) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp651 (NM_001166644) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Zfp651
Synonyms:	4732420M22Rik; BC030045; R74626
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG214800 representing NM_001166644
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTTGCTGGTTGAGAAGACGACAGACTCCCCGGCGGCAGAGTTCGCTGGTGGAGGACGTGGCCCTGC
ACTTCGCCTGCTTGATGGGCCGCCTCAATGAGCAGCGCCTCTCCAGCCAGACCTGTGCGATGTGGACCT
AGTGCTGGTGCCTCAGCATAGCGTGTTCGCCGCACACAAGGGTGTGCTGGCCGCCTACAGCCAGTCTTTC
CACTCGCTCTTACACAGAAACAAGCAGCTGCAGCGGGTCGAGCTATCTCTGGAGGCACTCGCGCCACGC
GCCTGCAGCAGATCCTCAACTTCATCTACACGTCCAAGCTGCTCGTGAACGCGGCCAACGTGCACAGGT
GCTCAGCGCGCCTCGCTGCTGCAGATGGCTGACATCGCGGCCTCCTGCCAGGAGCTGCTGGACGCTCGA
TCCTTGGCGCCCTCGGGCCCTGTGGCCCTGGCACAGCCTGCCACCAGCTGTGCTCCGGTGCCTCCACCAC
CCTACTACTGTGACATCAAGCAGGAGGCAGATGCCCCAGGCCTCCCAAAGATCTATGCCCGGAGGGGCC
CGATCCCTACTCTGTGCGTGTGGAGGATGGAGCAGGGGCTGCTGGGGATGTAGGCTCTGCCACTGCTGGG
CCAGCACAGACACTCTTCAAGGAGGAGAAGGAAGGGGCCCTGAAGAGGCAGAGGCCCCAGGTAGCCTGT
GCAAAGTGGAAAGTGGAGAGGGGCTGGAGCCAGAACTGGATGCCTCGGGACCTATGGGCACCAGGAGCA
GTCACAAATCATTGTTGAGGTGAACCTTAAACAACCAAACTGCATGTGTCCACCGGACCCGAGGGCAAG
CCAGGCTCTGGTGCCAAATCCAGCCACTGTGGTGTCTCGCCAGGAGGATGGCATGCAGGGACACTCGGAAG
AAGAGGAGGAGGAAGGAGGAGGGAGTGGAGGGGAGAGGAAGAGGAGGAAGAGGAAGAAGAAGAGGAAGG
CAGCCAGGAGGAGGAGGAGGAGGAAGAGGAGGAAGGGCCAGCGATCATCGAGAGGATGAGGATGAAGAA
GATGGGCCAGTGAGCAGGATGCAGAGAGTCTGAGGAGGAGAGGGAGGCTGAGGGCAGGCAGGACCCCG
CAGGTCCCGCTGGGTGTCAGGGCAGCCAGGTGGACCCTCCACCACATAGCCGAATGTCCACACGGTCCCG
GGGACAGAACACCCGGCCCGGGCCACCCTGAGCCAGAGGAGGCTGGACGTAGAGGTGGGAAGAGGCC
AAGGCCTCTGGAGCTGTTCCGGCATCCCAGGCAGCTGACGGGCTGGGGCCAAAGTGAAGCTGGAGGAAA
AACAGCAGCATCCGTGCCAGAAGTCCCTCGGGTCTTCAACAACCGCTGGTACCTAGAAAAGCATATGAA
CGTGACCCACAGCCGGATGCAGATCTGTGGCCAGTGTGGGAAGCGCTTCTGCTGGAGAGTGAGCTGCAG
CTGCACAGACAGACGGACTGTGAGCGCAACATCCAGTGCATGACATGTGGCAAAGCCTTCAAGAACTGT
GGTCTCTCATGAGCACAACAAGATAGTGCACGGCTACGCGGAGAAGAAGTCTCCTGTGAGATCTGCGA
GAAGAAGTTTACACCATGGCGCAGTGCACAGCAGATGGTTGCCACACCAAGGACATGCCCTTACC
TGTGAGACCTGTGGGAAGTCCTTCAAGCGCAGTATGTCGCTGAAGTCCACTCGCTGCAGCACTCGGGG
AGAAACCTTCCGCTGTGAGAACTGCAACGAACGCTTCCAGTACAAGTACCAGCTGCGGTACACATGAG
CATCCACATTGGCCACAAGCAGTTCATGTGTGTCAGTGGTGTGGCAAGGACTTCAACATGAAGCAGTACTT
GACGAGCAGATGAAGACACACAGGGGAGAAGCCTTATCTGTGAGATCTGTGGCAAGAGTTTACCA
GCCGCCCAATATGAAGCGGCATCGACGCACACACAGGAGAGAAGCCATACCCGTGTGACGTGTGCGG
CCAGCGTTTCCGCTTCTCCAACATGCTGAAGGCACACAAAGAAAAGTGTCTCCGTGTGAGCCACCCGCTG
CCCGGTGACCCCGCCACCCTGCCTGCCACGCATCTGCAGCCACAGCTCCGCTCTCCCCACAGCACCTC
CCAGGCTGGACACTAAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG214800 representing NM_001166644
 Red=Cloning site Green=Tags(s)

```

  MLLVEKTTDSPAAEFSLVEDVALHFACLMGRLENEQRLFQPDLCDVDLVLVPQHSVFPAAHKGVLAAYSQFF
  HSLFTQNKQLQRVELSLEALAPSGLQQILNFIYTSKLLVNAANVHEVL SAASLLQMADIAASCQELLDAR
  SLAPSGPVALAQPATSCAPVPPPPYYCDIKQEADAPGLPKIYAREGPDYSVRVEDGAGAAGDVGSATAG
  PAQTLFKEEKEGAPEEAEAPGSLCKLESGEGLEPELDASGTYGHQEQSQIIVEVNLNNQTLHVSTGPEGK
  PGSGANPATVVVLGQEDGMQGHSEEEEEEGGSGGGEEEEEEEEEGSQEEEEEEEEEGPSDHREDEDEE
  DGPSEQDAESSEEREAEGRQDPAGPAGCQGSQVDPPPHSRMSTRSRGNTRRRATPEPEEAGRRGKRP
  KASGAVPASQAADGLGAKVKLEEKQQHPCQKCPRVFNRRWYLEKHMNVTHSRMQICGQCGRFLLESELQ
  LHRQTD CERNIQCMTCGKAFKWLWLSHEHNKIVHGYAEKKFSCEICEKKFHTMAHVRKHMVAHTKMPFT
  CETCGSFKRSMMLKVHSLQHSGEKPFRCENCNERFQYKYQLRSHMSIHGHKQFMCQWCGKDFNMKQYF
  DEHMKTHTGKPYICEICGKSFTSRPNMKRHRRTHTGKPYPCDVCGRFRFSNMLKAHKEKCFRVSHPL
  PGDPATLPATHLQPTAPLFPTAPPRLDTN
  
```


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:


ACCN: NM_001166644

ORF Size: 2187 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166644.1](#), [NP_001160116.1](#)

RefSeq Size: 6353 bp

RefSeq ORF: 2190 bp

Locus ID: 270210

Cytogenetics: 9 F4