

Product datasheet for **MG214958**

1700024G13Rik (NM_001034037) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700024G13Rik (NM_001034037) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: 1700024G13Rik
Synonyms: C10orf53
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG214958 representing NM_001034037
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCTGCGCAAGCGGTGGTGACCTTGCCTATGGGCCCTACAGCGCAGTGGCCTGTCCGTGGAGCACC
 GCACCTACCGACTGCAGGGTCTGCAAGCTGTGTTGGCCAAAGACGGACATCAGATCATCTGGAGCAGAT
 AGAAGACTGGAATCTGGTGGAGCTTGTGGTAAATGAGGAGACGGTCTTCAATGTGATATCCAAGAGCTG
 GAGTTTGGGGTGATGGCAAAC TAGACCCACTGTGTAAGAGGCCAGGATCGCTGTGCTAAATGCCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG214958 representing NM_001034037
 Red=Cloning site Green=Tags(s)
 MPAQAVVTLRYGPYSAVGLSVEHRTYRLQGLQAVLAKDGHQIILEQIEDWNLVELVVNEETVFCQDIQEL
 EFGGDGKLDPLCEEARIAVLNAF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001034037

ORF Size: 279 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001034037.1</u> , <u>NP_001029209.1</u>
RefSeq Size:	399 bp
RefSeq ORF:	282 bp
Locus ID:	67085
UniProt ID:	<u>Q3KNL4</u>
Cytogenetics:	14 B