

Product datasheet for **MG215027**

Tspan8 (NM_001168680) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Tspan8
Synonyms:	C76990; E330007O2IRik; Tm4sf3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG215027 representing NM_001168680
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCAGGTGTGAGCAGCTGTTAAAGTATTCTATGTTTTCTTCAACTCTTGTCTGGGTATGTGGTA
CACTGATCTTGGGACTGGCCATATGGGTGAGAGTAAGCAAAGATGGAAAAGAGATCATCACGTCAGGAGA
CTCTAGCACGAACCCCTCATTGCTGTGAATATCTTGATTGCAGTTGGGTCTATCATCATGGTTCTAGGT
TTCCTGGGATGCTGTGGAGCTGTGAAAGAAAGCCGGTGCATGCTTCTCTGTTTTTTATCGGCTTGCTGC
TGATTCTGATTCTGCAAGTGGCAGCAGGCATCCTAGGAGCCGCTTCAAACCTGAGTATAACCGAATTCT
GAATGAAACGCTATATGAAAATGCTAAGCTTCTGTCGGACAACACTGATGAGGCAAAAAGACTCCAGAAG
GCCATGATTGTGTTTCAGTCGGAGTTCAAGTGTGGCTTGGAGAATGGAGCTGCCGACTGGGAAATA
ATTTTGTGGAGGCAAAAAGAGTCATGTCAATGTACAGGAAGTACTGTGCAACTATCAGGGGAGTCCGT
TTACCCAAAGACTGTCTTCTTTGATCAAAGATCTGTTGAAAAGAACATTATTATTGTCATTGGAATT
GCTTTTGGACTGGCCGTTATTGAGATTCTCGGTTTGGTGTCTTCTATGGTGCTATACTGCCAGATTGGAA
GCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >MG215027 representing NM_001168680
 Red=Cloning site Green=Tags(s)

MAGVSSCLKYMSMFFNFLFWCGTLILGLAIWVRVSKDGKEIITSGDSSTNPFIAVNILIAVGSIIIMVLG
 FLGCCGAVKESRCMLLLFFIGLLLILILQVAAGILGAAFKPEYNRILNETLYENAKLLSDNTDEAKDFQK
 AMIVFQSEFKCCGLENGAADWGNNFVEAKESCQCTGTDCA TYQGSSVYPKTCLSLIKDLFEKNIIIVIGI
 AFGLAVIEILGLVFSMVLVCQIGSK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001168680

ORF Size: 705 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001168680.1](#), [NP_001162151.1](#)

RefSeq Size: 1416 bp


RefSeq ORF: 708 bp

Locus ID: 216350

UniProt ID: [Q8R3G9](#)

Cytogenetics: 10 D2

Product images:



Circular map for MG215027