

## Product datasheet for **MG215086**

### Olf1015 (NM\_146571) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Tag:** TurboGFP  
**Symbol:** Olf1015  
**Synonyms:** MOR213-3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG215086 representing NM\_146571  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGGATGTAGACAATCGCACCATCCTGACTGAATTTATTCTTCTGGGCTTCTCAGCAGACCCCATGGC  
 AGCTGATTCTATTTGGAATATTTCTAACCATCTACTTGATGACCTGTTAGGGAACATGACCTTGATTAT  
 TTTAATCCGCATCGATTCCAGGCTGCACACCCATGACTTTTTTATTGGAGGTCTATCTTTTTGGAC  
 TTCTGGTATAATTCTGTATACATACCAAAAATCTGGTCAACTGTGTTTCAGAGGATAAGCGTATTTCTT  
 TGGCTGGTTGTGGGCTCAGTTTTTCTTTTCTGTGTAGCAGCCTACACTGAGTGCTACCTGCTGGCAGC  
 CATGGCATATGATCGACATGCTGCAATTTGTAGCCATTGCTCTATTCAAGCATCATGTCCACTTCTCTC  
 TGTGCTGGGCTAGTGGGGCCTCCTATGTGGGAGGTTTTTGAATGCCATAGCCATACTGCCAACACTT  
 TCAGGCTAAGGTTCTGTGGTAAAAATATTATTGATCACTTTTTCTGTGATGTGCTTCCATTGGTAAAAAT  
 GTCTTGACAGACCCCGTGTCTATGTGAAAATCCTTTCCAGTATGGTGGGCTTCACTGTCCCTCAAGC  
 ATTCTTGCCATCCTAATTTCTATCTCAACATCCTGCTGGCTATACTGAGGATCCGCTCAGCCTCAGGAA  
 GACGAAAGGCTTTCTCCACCTGTGCCTCTCACCTGGTCTCTGTTACACTCTTCTATGGCTCCCTGCTCTT  
 CATGTATTCAAGACCCAGTTCCAACACTCCCTGGAGAGAGACAAAGTAGCTGCCATGTTCTATATCATC  
 ATCAATCCATTAACAACCTTTTATCTATAGTTTGAGAAAACAAAGATGTCAAAGAAGCCTTTAAGAAGT  
 TGATGCAGACTATAAAACAGCAAACC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



**Protein Sequence:** >MG215086 representing NM\_146571  
Red=Cloning site Green=Tags(s)

MDVDNRTILTEFILLGFSADPHWQLILFGIFLTIYLMTLGNMTLIIILIRIDSRLHTPMYFFIGGLSFLD  
 FWYNSVYIPKILVNCVSEDKRISLAGCGAQFFSCVAAYTECYLLAAMAYDRHAAICSPLLYSSIMSTSL  
 CAGLVGASYVGGFLNAIAHTANTFRLRFCGKNIIDHFFCDVLPVKMSCTDTRVYVKILSSMVGFTVLSS  
 ILAILISYLNILLAILIRISASGRKAFSTCASHLVSVTLFVGSLLFMYSRPPSSNYSLERDKVAAMFYII  
 INPLLNPFIYSLRNKDVKEAFKMLMQTIKQQT

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**ACCN:** NM\_146571

**ORF Size:** 936 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<a href="#">NM_146571.1</a>
<b>RefSeq Size:</b>	1016 bp
<b>RefSeq ORF:</b>	939 bp
<b>Locus ID:</b>	258564
<b>Cytogenetics:</b>	2 D
<b>Gene Summary:</b>	<p>Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]</p>

Product images:



Circular map for MG215086