

Product datasheet for **MG215153**

Slc24a2 (NM_172426) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc24a2 (NM_172426) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Slc24a2
Synonyms:	2810021B17Rik; 6330417K15Rik; AI847460; Nckx2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MG215153 representing NM_172426
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATCTCCACCAAAGCCCCACAGCCAGGCTCCTGCAGAAATGGTGTTCCCATGAGTCACCATTTGGCT
 GCAGGAGACATTATAATAGCAGGAAAAAACTGAAGCTAATTCGAGTCATTGGCCTTGTCTATGGGCCTGGT
 AGCTGTGAGCACTGTCCATTTTCAATCAGTGCCTTTACTGAGACAGACTCTCAGAGCAACAGAGGAGAA
 GCCAGTGATATGAGTGGCCCCAGAGTAGCACAGGGTCACCGTCAAAGAACTCTCTTAGACTTAAATGACA
 AGATTCGAGATTACCCCCACAGCCACCTGCTTCTCAGGAAGACCAAGCTGAGAACAGTACTGAGCACAC
 CCAGGGGAGACTACCCAAAGGACATCTTTTCCCTGGAGGAGCGAAGGAAAGGTGCCATCATTCTGCATGTC
 ATTGGGATGATCTATATGTTTCATAGCATTAGCCATCGTCTGTGATGAGTCTTTGTTCTCTCTAACTG
 TCATCACTGAAAACTGGGAATCTCTGATGATGTGGCTGGAGCCACCTTCATGGCTGCAGGCGGGTCAGC
 CCCAGAACTTTTACGTCTCTTATAGGGGTCTTCATTGCTCACAGCAATGTCGGCATTGGCACCATCGTG
 GGTTCCGCCGTGTTCAATATTCTCTTTGTTATTGGCATGTGTGCTCTGTTTTCTAGAGAAATCCTAAACC
 TGACCTGGTGGCCGCTCTTTCGGGATGTGCTCTTCTACATCGTTGATTGCTCATGTTGATCACATTTTT
 CCTGGATAACGTTATCATGTGGTGGGAAAGCTTACTCCTCTTAACAGCTTATTTGCGCTATGTGGTTTTT
 ATGAAGTTCAACGTCCAAGTAGAGAGATGGGTGAAGCAATGATATCTCGCAATAATGTTATCAAGGTGA
 CAGTGCCAGAAGCCCAAGCAAAGTACCTACAGCTGGGGACAAGGATGGGCCAACTTTACCGTCGAAGCC
 ACGCCTCCAACGAGGCGGAAGCTCTGCCTCCCTCCACAACAGTCTCATGAGGAATAGCATCTTCAGCTC
 ATGATACACACCTTGACCCACTCGCTGAAGAATTGGATCATATGAAAACTAAATATTATGACACAA
 TGACTGAAGAAGGGAGGTTTCAGAGAAAAAGCTTCAATTCTCCACAAGATCGCCAAGAAGAAATGCCAAGT
 GGATGAGAATGAGAGGCAGAATGGGGCTGCCAATCATGTGAAAAAATTGAACTCCCAACAGCACAAAGC
 ACAGAGGTGGAATGACGCCATCGAGCGAAGCTTCGGAACCCGTCGAGAATGGAACCTCTCCACAACA
 TTGAAGCTGCAGATGCCCGAAGGCCACAGAGACTGCCGAGGAAGAAGATGACCAGCCCCCTCAGCTTGTC
 CTGGCCGACCAACACCCGCAAGCAAGCCACCTTCTAATCGTTTTTCCATTGTGTTCCCTCTCTGGATC
 ACCTTACCTGACGTCCGAAAACCTGCTTCCAGGAAGTTCTTCCCTATCACATTCTTCGGCTCCATCACCT
 GGATTGCAGTTTTCTCATACCTGATGGTGTGGTGGGCACACCAGGTTGGAGAGACAATCGGCATTAGTGA
 AGAAATCATGGGTCTGACCATCTTGGCTGCTGGGACCTCTATCCCTGATCTTATCACCAGTGTATAGTG
 GCCCGAAAGGGCCTGGGGACATGGCTGTGTCCAGCTCTGTTGGAAGCAACATCTTTGACATCACTGTGG
 GGCTCCCACTGCCCTGGCTCCTCTACACCATCATCCACAGATTACGCCAGTGACTGTGACGAGCAATGG
 CCTCTTCTGTGCCATCGTCCTTCTCTTCATCATGTTGCTTTTGTCAATTCTCTCCATCGCCCTCTGCAAA
 TGGCGGATGAACAAAATCCTCGGCTTATCATGTTTGGCCTCTACTTCGTGTTCTGGTGGTCAGCGTCC
 TCCTGGAAGACAAAGTTCTCGTGTGCCCTGTCTCCATC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG215153 representing NM_172426
 Red=Cloning site Green=Tags(s)


MDLHQSPARLLQKWCSESHPFGCRRHNSRKKLKLIRVIGLVMGLVAVSTVPFISISAFETETDSQSNRGE
 ASDMSGPRVAQGHQRRTLLDLNDKIRDYTPQPPASQEDQAENSTEHTQGDYPKDIFSLEERRKGAILHV
 IGMIYMFIALAIVCDEFFVPSLTVITEKLGISDDVAGATFMAAGGSAPELFTSLIGVFIHNSVNGIGTIV
 GSAVFNILFVIGMCALFSREILNLTWWPLFRDVSFYIVDLLMLITFFLDNVIMWWESLLLLTAYFAYVVF
 MKFNVQVERWVKQMSRNNVIKVTYPEAQAKSPTAGDKDGPLPSKPRLQRGGSSASLHNSLMRNSIFQL
 MIHTLDPLAEELGSYGKLYYDTMTEEGRFREKASILHAKKKQVDENERQNGAANHVEKIELPNSTS
 TEVEMTPSSEASEPVQGNLSHNIEAADAPKATETAEEDDQPLSLSWPTNTRKQATFLIVFPVPLWI
 TLPDVRKPAKRFFPITFFGSITWIAVFSYLMVWWAHQVGETIGISEEIMGLTILAAGTSIPDLITSVIV
 ARKGLGDMAVSSSVGSNIFDITVGLPLPWLLYTIHRFSPVTSSNGLFCAIVLLFIMLLFVILSIALCK
 WRMNKILGFIMFGLYFVFLVSVLLEDKVLVCPVSI

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_172426

ORF Size: 1998 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_172426.2](#), [NP_766014.1](#)


RefSeq Size: 10606 bp

RefSeq ORF: 2001 bp

Locus ID: 76376

Cytogenetics: 4 C4

Product images:



Circular map for MG215153