

## Product datasheet for **MG215163**

### **Pira1 (NM\_011087) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pira1 (NM_011087) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pira1
Synonyms:	6M21; Ly89; Pir; PIR-A1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MG215163 representing NM\_011087  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCCTGCACCTTCACAGCCCTGCTCCGTCTTGACTGACTCTGAGCCTCTGGATCCCAGTGCTGACAG  
GGTCCCTCCCTAAGCCTATCCTCAGAGTACAGCCAGACTCTGTGGTCTCCAGGTGGACTAAGGTGACTTT  
CTTTTGTGAGGAGACAATTGGAGCCAATGAGTACCGCCTCTATAAAGATGGAAAGCTATATAAACTGTA  
ACAAAGAACAACAGAAGCCAGCAAACAAGGCTGAATTCTCACTCTCAAATGTAGACCTGAGAAATGCAG  
GTCAATATCGATGTTCTACAGCACCCAGTATAAATCATCAGGCTACAGTGACCCCTGGAGCTGGTGGT  
GACAGGAGACTACTGGACACCCAGCCTTTAGCCCAAGCCAGCCCTGTGGTAACTTCAGGAGGGTATGTC  
ACCCTCCAGTGTGAGTCTGGCACAACGATCACAAGTTTATTCTGACTGTAGAAGGACCACAGAAGCTCT  
CGTGGACACAAGACTCACAGTATAATTACTCTACAAGGAAGTACCACGCCCTGTTCTCTGTGGGCCAGT  
GACCCCAACAGAGATGGATATGCAGATGTTACAGTTATGACAGGAACAGACCATATGTGTGGTCACT  
CCAAGTGAATCCGTGGAGCTCCTGGTCTCAGGTAATCTCCAAAAACCAACCATCAAGGCTGAACCAGGTC  
CTGTGATCGCCTCCAAAAGAGCAATGACCATCTGGTGTCAAGGGAACCTGGATGCAGAAGTATATTTTCT  
GCATAATGAGGGAAGCCAAAAACACAAAGCACACAGACCCTACAGCAGCCTGGGAAACAAGGCAAGTTC  
TTCATCCCTTCTATGACAAGACAACATGCAGGGCAATATCGCTGTTATTGTTACGGCTCAGCTGGTTGGT  
CACAGCCAGTGACACCTGGAGCTGGTGGTGACAGGAATCTATGAACACTATAAACCCAGGCTGTCAGT  
ACTGCCAGCCCTGTGGTGACAGCAGGAGGAAACATGACACTCCACTGTGCCTCAGACTTCACTACGAT  
AAATTCATTCTCACAAGGAAGATAAGAAATTCGGCAACTCACTGGACACAGAGCATATATCTTCTAGTA  
GACAGTACCGAGCCCTGTTTATTATAGGACCCACAACCCCAACCCATACAGGGACATTGATGTTATGG  
TTACTTCAAGAATGCCCCACAGCTGTGGTCACTACCTAGTGATCTCCAACAAATACTCATCTCAGGGCTG  
TCCAAGAAGCCCTCTCTGCTGACTACCAAGGCCATATCCTGGACCTGGAATGACCCTCACCTGCAGT  
GTTTCTCTGACATCAACTATGACAGATTTGCTCTGCACAAGGTGGGGGAGCTGACATCATGCAGCACTC  
TAGCCAGCAGACTGACACTGGCTTCTCTGTGGCAACTTCACACTGGGCTATGTGAGTAGCTCCACTGGA  
GGCAATACAGATGCTATGGAGCACACAACCTTCTCTGAGTGGTCAGCCTCCAGTGGCCCTGGACA  
TCCTGATCACAGGACAGTCCCTCTCACTCCTTCCCTCTCAGTGAAGCCTAACACACAGTGCAGTCAAG  
AGAGACCGTGAGCCTGCTGTGTTGGTCAATGGACTCTCTGGATACTTTCATTCTGTCCAAGGAGGGATCA  
ACCCAGCAACCCCTACAATAAAATCAAAGTCCCATGATCAGCAGTCCCAGGCAGAATTCCATGAGTG  
CTGTGACCTCCCATCTCTCAGGCACCTACAGGTGATGGTGCTCAAGACTCATCTTTCTACCTCTGTGC  
ATCTGCCAGTGCCCTGTGGAGCTCACAGTCTCAGAAACCATTGAATCCTCCACCTGGTCTCCAAAAGG  
CCCATCCCAATTCCACAGAGAACAAGGATCACACAATGGAGAATCTCATCAGGATGGGGATGGCAG  
TCCTGGTCTCATCGTCTTTGATTCTAGCCACTGAGGCTTGGCGAAGCCATAGACAGACCCACCTGC  
AGCTGGGAAC

**ACGCGT**ACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG215163 representing NM\_011087  
 Red=Cloning site Green=Tags(s)


MSCTFTALLRLGLTSLWIPVL TGLPKPILRVQPDSSVSRWTKVTFCEETIGANEYRL YKDGKLYKTV  
 TKNKQKPANKAEFSL SNVDLRNAGQYRCSYSTQYKSSGYS DPLELVVTGDYWTPSLLAQA SPVVTSGGYV  
 TLQCESWHNDHKFILTVEGPQKLSWTQDSQYNYSTRKYHALF SVGPVTPNQRWICRCYSYDRNRPYVWSP  
 PSESVELLVSGNLQKPTIKAEPGPVIA SKRAMTIWCQGNLDAEVYFLHNEGSQKTQSTQTLQQPGNKGF  
 FIPSMTRQHAGQYRCYCYGSAGWSQPSDTLELVVTGIYEHYKPRLSVLPSPVVTAGGNM TLHCASDFHYD  
 KFILTKEDEKFGNSLDTEHISSSRQYRALFIIGPTTPTHTGTFRCYGYFKNAPQLWSVPSDLQQIILISGL  
 SKKPSLLTHQGHILDPGMTLTLQCFSDIN YDRFALHKVGGADIMQHSSQQTDTGFSVANFTLGYVSSSTG  
 GQYRCYGAHNLSSEWSASSEPLDILITGQLPLTPSL SVKPNHTVHSGETVSLLCWSMDSLDTFILSKEGS  
 TQQPLQLKSKSHDQQSQA EFSMSAVTSHLSGTYRSYGAQDSSFYLLSSASAPVELTVSETIESSTWSPKR  
 PIPPIPTENKDHTMENLIRMGMAVLVLI VLSILATEAWRSHRQTHPAAGN


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



**Plasmid Map:**


**ACCN:** NM\_011087

**ORF Size:** 2040 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_011087.1](#), [NP\\_035217.1](#)

**RefSeq Size:** 3453 bp

RefSeq ORF: 2043 bp  
Locus ID: 18722  
Cytogenetics: 7 2.15 cM