

Product datasheet for **MG215168**

Zfp597 (NM_001033159) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Zfp597
Synonyms:	4933407K12Rik; AI462383; AI666357
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG215168 representing NM_001033159
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGTCCACGCTCCCCACATCGGATGAGCAGAGACCACTGCTATTTGAAGATCTGGATGTGTATTTTT
CTCAAGAGGAATGTGTGAGTCTGCACCCTGCTCAGAAGACCCTCAGCAGAGAAGCCCCATTGAAGAGTTT
TGAGAAATTTGGACTCGATAGGAGGGGAAGACCAGACTGAAAGCAGTCAGGAGTTGAGCCTCGAGCCTGTG
GAACCTGAAGAGCTGTCTCTAGAAAAGGACTCCATTGCTGCACCCGTTGTTCACTACTCTGAGCAATCTG
AGCAGGATGATGAAATACTTGAAGACAAAATGTCAGGTGGAACGTTGACTTGCAAGACAAACTTATAAG
CCTCTTGTTACCATTTGACAACCAACCCCATAGTAGAATTATCTCAATGTTTAGGAGTAAAAACACTT
TCAGATATTATTGAAGTTCCTTGGGAAGAAGCCAAAAATGTGTACAAGTGTCTGACTGTGACCAAAAGCT
TCAGTGATAATACATACCTTGTTCATCAGAAAATTCATTTAAGAGAGAAAAAGTATAAATGTAGTAC
CTGTGAGAAGACCTTCAGTCACCGAACCAACCTGAGGACACATAAGCGAATCCACACTGGTGAAGAGCCT
TACAAGTGTACCGAGTGTGCTGCCAGCTTCCGACAGCAGTACACCTGTCTAGGCACATGAATAGCCATT
TAAAGGAGAAAATATACATGTAGTGTATGTGGGAGAGGTTTTATGTGGCTCCAGGATTGGCAGAACA
TCAGAAGAGTCATACTGATAAGAAATCTTATGAATGTGCTGACCACAATCGGAAGCAAATCTGGCTTTG
CCTGAAGAGAGAGGTTTCATCAGACACACTGCCCCAGCACACACACTGTGTGAAGACCTTGGAGCAGCCCT
CAGACCCTACTCTCCCTGAGAAAGACCATAAGGAGGACTCGAAAACTGCAGTATTGATGATGAGGACTT
TTTTTCATTCTCCAGATTCAAACCTTACAGTGTCTTGATTGTGACATGACTTTTTCTTGTTCTCGGAG
CTTGTCTCTACCAAACCATTCATGACCTGAAAAACCTCATAAGTGCAAAACATGTACAAAACCTTTTG
CTTTGAGTCAGAGCTTGTATCCCATGAGAAGGCCACAGAAGAGAAGGCCTTTTAAATGCACAGTGTG
TGGAAGAGTTTCAGTGTGAGTATGCGACTCGTCACTCATAACGAGCCATAGGAGAAACAGCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >MG215168 representing NM_001033159
Red=Cloning site Green=Tags(s)

MASTLPTSDEQRPLLFEDLDVYFSQEECVSLHPAQKTL SREAPLKSFENLDSIGGEDQTESSQEL SLEPV
 ELEELSLEKDSIAAPVVHYSEQSEQDDEILEDKMSGGTLTKTKLISLLVTIDNQTPLVLSQCLGVKTL
 SDIEVPWEEAKNYYKCPDCQSFSDNTYLVLHQIHLREKKYKCSCEKTF SHRTNLRTHKRIHTGEKP
 YKCTECAASFRQQSHLSRHMNSHLKEKLYTCSVCGRGMWLPGLAEHQKSHTDKKS YECADHNREANLAL
 PEERGSDDLTPQHHTCVKLTLEQPSDPTLPEKDHKEDSKNCSIDDEDFFSF SRFKPLQCLDCDMTPFCFSE
 LVSHQTIHDLKPHKCKTCTKTF AFESELVSHEKSHRREEPFKCTVCGKSFVSMQLVTHKRAHRNSK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001033159

ORF Size: 1257 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001033159.2](#), [NP_001028331.2](#)


RefSeq Size: 5630 bp

RefSeq ORF: 1260 bp

Locus ID: 71063

Cytogenetics: 16 A1

Product images:



Circular map for MG215168