

Product datasheet for **MG215286**

Pbsn (NM_017471) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pbsn (NM_017471) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Pbsn
Synonyms: PB; Prbs
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG215286 representing NM_017471
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGATGAGGGTCATCATCTCTGCTCACATGCATGTGCTAGGCGTCTCCAGTGTGATGAGTCTCAAAA
AGAAGATTGACGGGCTTGCAAACAATTTACTTAGCTGCCAGTACCATGGAAAAGATAAATGAAGGCTC
ACCATTGAGAACCTACTCCGTCACATTTGTGTGGGAGGAGATGCAACCAAGTCTACCTCTATTTTTTT
ATTAAGAAAGGGACTAAGTGCCAACTGTATAAAGTCATAGGAAGGAAAAACAAGAAGTTTACTATGCAC
AGTATGAAGGGAGCATAGCATTATGTTAAAGATGGTAAATGAGAAGATATTGTTGTTTCATTATTTTAA
CAAGAACAGAAGAAATGATGTACACGAGTGGCTGGAGTTTTGGCAAAAGGCAAACTGAATAAGGAG
GAGATGACGGAGTTCATGAACCTAGTGAAGAAATGGGCATTGAGGAAGAGAATGTACAGCGTATCATGG
ACACGGACAACCTGTCCAAGCAAGATCAAACCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG215286 representing NM_017471
Red=Cloning site Green=Tags(s)

MMRVIIILLTLHVLGVSSVMSLKKKIDGPWQTIYLAASTMEKINEGSPLRTYFRHILCGRRCNQVYL YFF
 IKKGTKCQLYKVI GRKKQEVYYAQYEGSIAFMLKMVNEKILLFHYFNKNRRNDVTRVAGVLAKGQLNKE
 EMTEFMNLVEEMGIEEENVQRIMDTDNCP SKIKP

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



View online »

Cloning Scheme:



ACCN: NM_017471

ORF Size: 522 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_017471.2](#), [NP_059499.1](#)

RefSeq Size: 805 bp


RefSeq ORF: 525 bp

Locus ID: 54192

UniProt ID: O08976

Cytogenetics: X B

Product images:



Circular map for MG215286