

Product datasheet for **MG215300**

Gm4836 (NM_009529) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: Gm4836
Synonyms: 22526
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG215300 representing NM_009529
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTCTATTAAGAACTGTGGGTGATACCAAAGGATGGTTACTTATTACTTCTTGACTATGATTCAGATG
 AAGAAGAAGAGCAGGCTCATTCTGAAGTAAAGAGGCCAGCATTGGCAAACATGAGAATATGCCGCCTCA
 CGTAGAAGCAGATGAAGATATAAGAGATGAACAAGACAGTATGTTGGATAAACTGGAGAAAACGTAAGT
 TTCTCAGAGGAATGGCAGCGTTTTGCACGTTCTGTAGAGACACCAATGGAAAACGGAATTTGTTAAGTG
 GTGAACAGCAAGTTAGGAATGCTTCAGAATTGGACCTTATGGAAGTACAGAATCCAGTAACTCATGATGA
 TGGGAATGCAAATCCTGAAGAAGTAGTTGGAGATACACGAAAGAAGATCAACAACAAATTTGTGAGCAG
 AAGTTTGATATGGATACAGAAATCAATGAAGAACAAGAAAAATCAGTGAACAATTATCAAAAAGAAC
 AACAAAGCATTAAAAGTCCGAATGTAGTCAGAGTCCGACCATGGAAGCAATTGAAGACATGCATGAGAA
 GTCCATGGAGGTTTGATGAACATGGAGACCAACAACACTACGATATGCTTTTTGATGTAGATGGTGAAGAG
 ACATTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG215300 representing NM_009529
 Red=Cloning site Green=Tags(s)


MSIKKLWVIPKDGYL LLLDYDSDEEEEQAHSEVKRPAFGKHENMPHVEADEDIRDEQDSMLDKSGENV
 FSEEWQRFARSVETPMENWNLLSGEQQVRNASELDLMEVQNPVTHDDGNANPEEVVGDTRKKINNKLC
 EKFDMDIQKFNEEQEKSVNNYQKEQQALKLSECSQSPTMEAIEDMHEKSMEGLMNMETNNYDMLFDV
 DGEE
 TL

TRTRPLE - GFP Tag - V



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_009529

ORF Size: 636 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_009529.3](#), [NP_033555.1](#)

RefSeq Size: 917 bp
RefSeq ORF: 639 bp
Locus ID: 22526
Cytogenetics: X

Product images:

Circular map for MG215300