

Product datasheet for **MG215326**

Kir3dl2 (NM_177748) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kir3dl2 (NM_177748) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Kir3dl2
Synonyms:	Kirl2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG215326 representing NM_177748 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGCTCTGGTTCCTCAGCCTCATGTGTTCTGGTTTCTCCTAGTCCAGAGGATGTCAGCACACGTGG
GAAGTCATGATAAACCATTCTATATGCCTGGCCAAGTTATGTTGTTCTCTGGGACAGAATGTGACTCT
TACATGTGACTCTCATAGAGGATCTAACATATTCAAACTGTACAAAGAAGAAGGGTCTCCAATCCCTCAG
TTACATGAGACTACATTCAGAAGAGCCAGGTCTTTGGTCTGTGACCACAGAACATGCAGGAACCTACA
GATGCTCCATCCTCAATATGCTAATGTATTGTCAGCACATAGTGAGCCTCTGAAAATCATAATTTCAGG
AATCTACTTGAAGCCTTTCCTCTTGATTCTACAAAGTCCCTGGTGAATTCAGGAGGGAATGTGACCCTG
GAATGTCACCTCAGAGAATATGTTTCGATACATACATTTTGATTTCCACAGAAATGGGAATAATCAAGAACT
CTGTCCAAGTTTCTGCAGAACATCATGAAAGTGGGTACATGTCACCTACTCAATAGGTCCCATGACACC
TGACCTTGTGGGACTTACACATGCTATGGTGCTAACAGTTACTACCCCTATGAATGGTCAGACCCCACT
GACCCATTGACATCAAGATCACAGGTGTATACAAGAAACCTTCTCTGTGAGTCTGTATGGGGCCTGTGC
TGATGATGTCAGGAGAAACCATGACCCTGTCTTGATCTCTGACCATCAGTTTGACATGTTCCATATGTC
CAGGGAAGGGTGCCTCAGGGACAGGGGATGCCTGGAGTGCAGATCCACAGTGGTAAATTCGAGGCAAAA
TTCCTTTAGTTCAATGATCCAAAAGGGGAACATAGATGCTATGGCTTTTCAGGAACCTTCCCATG
TGTGGTCTCCTCCCAAGTACCCCTGTACCTTCTGCCAAAGGAACTGTCTGCATGCACAGAAGAAGA
TCCAAAATACATAACTGCAAGAACCTGCGCATTCTAACTGGGCTCTTAATGACCATGGTCTTGTGCATC
ATCATCATCTTTTATTCTTGTACTTTCCCAAACAGAATAAGTCTCAGAAACAAGCAGCAGCCAGCATGG
AACAAAGAGTCTGAAGTTAAAAACACAATAAACAGACAGCTTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG215326 representing NM_177748
 Red=Cloning site Green=Tags(s)

MLLWFLSLMCSGFFLVQRMSAHVGSHPKFLYAWPSYVVPLGQNVTLTCDSHRGSNIFKLYKEEGSPIQ
 LHETTFQKSQVFGPVTTEHAGTYRCFHPQYANVLSAHSEPLKIIISGIYLPFLILQSPLVNSGGNVTL
 ECHSENMFDTYILISHRMGIKNSVQVSAEHHESGSHVYSIGPMPDLVGTYTTCYGANSYYPYEWSDPS
 DPIDIKITGVYKPPSLVLMGPVLMMSGETMTLSCI SDHQDFMFHMSREGVPQGQMPGVQIHSKGFEAK
 FLLSSMIQKGNRYCYGSRNSSHVWSSPDPLYLPAKGNCPACTEEDPKIHNCNLRILTGLLMTMVLVI
 IIFYSYCFPKQNKSQKQAAASMEQESEVKKHNKQTA

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja1823_e10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_177748

ORF Size: 1164 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_177748.2, NP_808416.1</u>
RefSeq Size:	1734 bp
RefSeq ORF:	1167 bp
Locus ID:	245615
Cytogenetics:	X 58.45 cM

Product images:


Circular map for MG215326