

Product datasheet for **MG215339**

Rtf1 (NM_030112) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rtf1 (NM_030112) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rtf1
Synonyms:	2900005O08Rik; 6530416A09Rik; AI853581; AU043144; AW553985; cbp82; Gtl7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG215339 representing NM_030112
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCGCGGTGCCTTTGTGTGGTTCGAGCGCGCGCGGTGGCGCGCGCGGTGGCGCGCGCGCAGTGGCGG
 TCCCACTGGCAGGCGGGCAAGAGGGGAGTCAGGGCGGCGTCCGGCGTGAAGCCGGGGGACCACCATGGT
 AAAGAAGCGAAAAGGCCGTGTCGTGATCGACTCGGATACGGAAGACAGCGGCAGTGACGAGAACCTGGAT
 CAGGAGCTCTTGTCCCTGGCGAAAACGGAAGCGCAGTACTCCGAGGAGAAGGAGCCTCCTGTGAGTCAGC
 CTGCGACCTCATCAGATTCAGAGACCTCCGACAGTATGATGAGTGGACATTTGGGAGCAATAAAAAATA
 GAAGAAAGGAAAAACCAGAAAAGTAGAGAAGAAAGGAGCCATGAAGAAGCAAGCCAAACAAAGCTGCCTCC
 TCGGGCAGCTCAGACAGAGACAGCTCAGCTGAGAGCTCTGCCCTGAGGAAGGTGAAGTGCAGATTCTG
 AGAGCAGCAGCTCCTCTCCAGTTCAGATTCACTCCTCGTCAGAAGATGAAGAATTCATGATGGCTA
 TGGAGAGGATCTTATGGGAGATGAGGAAGATAGGGCCCGACTGGAGCAGATGACGGAGAAGGAGAGGGAG
 CAGGAGCTCTTCAACCGCATAGAGAAGAGAGAGGTGTTGAAAAGAAGATTTGAAATTAAGAAAAACTAA
 AAACAGCCAAAAAGAAAGAGAAGAAAAAGAAAAAGAAAAAGCAGGAAGAAGAAACAGAAAAAGAAAGCT
 GACACAGATTCAGGAATCTCAGGTTACTTCCCACAATAAAGAACGGCGTTCGAAGCGGGATGAAAACTA
 GATAAGAAATCTCAAGCCATGGAGGAGCTAAAAGCAGAGCGAGAAAAACGGAAGAACAGAACAGCTGAGC
 TCCTCGCCAAAAACAGCCATTGAAAACAGTGAAGTCTACTCTGACGATGAGGAGGAGGAAGATGATGA
 CAAGTCCAGTGAGAAATCGGACCGCTCATCCCGGACATCATCTGATGAAGAGGAAGAGAAAGAA
 ATCCCTCCAAAATCACAGCCAGTCTCATTACCCGAGGAAGTGAATCGGGTTCGGTATCACGGCATAAGC
 TAGAGCGCTGGTGTACATGCCCTCTTTGCTAAAAGTGTACAGGATGTTTTGTACGGATTGGCATTGG
 AAACCAACAATAGCAAGCCAGTTTACCGGGTGGCTGAGATCACGGGTGTTGTGAAAAGTCCAAAGTTTAC
 CAGCTGGGTGGCACCAGAACAAATAAAGGGCTCCAGCTACGGCATGGCAATGACCAGCGAGTGTCCGCC
 TAGAGTTTGTGTCCAACCAAGAATTCAGTAAAAGTGAATTCATGAAGTGGAAAGAAGCGATGTTTTCTGC
 TGGCATGCAGTTGCCACACTAGATGAAATCAATAAAAAGGAATTATCTATTAAGGAAGCTCAATTAT
 AAATTCATGATCAGGACATTGAGGAGATTGTAAGAGAGAAGGAAAGGTTCAAGAAAGCCCAACAACT
 ATGCTATGAAGAAGACTCAGCTGCTAAAGGAGAAGGCTATGGCTGAGGACCTGGGAGATCAGGACAAAGC
 CAAACAAATCAAAGACCAGCTAAATGAGCTAGAGGAGCGAGCAGAGGCCCTGGACCGCCAGCGGACCAAG
 AACATATCTGCTATCAGTACATCAACCAGCGAATCGGGAGTGAATATTGTGGAGTCTGAGAAGGCTC
 TGGTGGCTGAAAAGTCACAACATGAGAAACCAACAGATGGACCCCTTACCAGCGCGCAGTCAAACCCAC
 CATCGTCTCAATTCAGGGATCCAGCTGTTCCAGGCTGCCATTTTGGCCAGCTCAATGCAAAATACGGT
 TCTGGGTATTACCAGATGCTCCAAAGGAAATGAGCAAGGGTCAAGGAAAGGATAAAGATTTGAATTCTA
 AGACAGCCAGTGACCTCTCGAAGACCTGTTCAAAGTGCATGACTTTGACGTGAAGATTGACTTACAGGT
 TCCAAGCTCAGAGTCAAAGCTTTGGCCATTACCTCAAAGCCCCACCAGCCAAGGACGGAGCTCCGAGG
 AGATCTCTGAACCTAGAAGACTACAAAAACGACGGGGCTTATT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG215339 representing NM_030112
 Red=Cloning site Green=Tags(s)

```
MRGRLCVGRAAAVAAVAAAAVAVPLAGGQEGSQGGVRRGSRGTTMVKKRKGRVVIDSDTEDSGSDENLD
QELLSLAKRKRSDSEEKPPVSPQAASSDSETSDSDDEWTFGSNKNKKKGKTRKVEKKGAMKKQANKAAS
SGSSDRDSSAESSAPEEGEVSDSESSSSSSSDSDSSSEDEEFHDGYGEDLMGDEEDRARLEQMTEKERE
QELFNRIEKREVLKRRFEIKKKLKTAKKKEKKEKKEKKEKKEKKEKKEKKEKKEKKEKKEKKEKKEKKEK
DKKSQAMEELKAEREKRNRTAELLAKKQPLKTEVYSDDEEEEDDKSSEKSDRSSRTSSSDEEEKEE
IPPKSQPVSLPEELNRVRLSRHKLERWCHMPFFAKTVTGC FVRIGIGNHNSKPVYRVAEITGVVETAKVY
QLGGTRTNKGLQLRHGNDQRFVLEFVSNQEFTESEFMKWKEAMF SAGMQLPTLDEINKKELSIKEALNY
KFNDQDIEEIVKEKERFRKAPPNYAMKKTQLLKEKAMAEDLDGQDKAKQIQDQLNELEERAEALDRQRTK
NISAISYINQRNREWNIVESEKALVAESHMNRNQMDPFTRRQCKPTIVSNRDP AVQAAILAQLNAKYG
SGVLPDAPKEMSKGQGDKDLNSKTASDLSEDLFKVHDFDVKIDLQVPSSSKALAITSKAPPAKDGAPR
RSLNLEDYKRRGLI
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_030112

ORF Size: 2145 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030112.2](#), [NP_084388.2](#)

RefSeq Size: 4680 bp

RefSeq ORF: 2148 bp


Locus ID: 76246

UniProt ID: [A2AQ19](#)

Cytogenetics: 2 E5

Gene Summary: Component of the PAF1 complex (PAF1C) which has multiple functions during transcription by RNA polymerase II and is implicated in regulation of development and maintenance of embryonic stem cell pluripotency. PAF1C associates with RNA polymerase II through interaction with POLR2A CTD non-phosphorylated and 'Ser-2'- and 'Ser-5'-phosphorylated forms and is involved in transcriptional elongation, acting both independently and synergistically with TCEA1 and in cooperation with the DSIF complex and HTATSF1. PAF1C is required for transcription of Hox and Wnt target genes. PAF1C is involved in hematopoiesis and stimulates transcriptional activity of KMT2A/MLL1. PAF1C is involved in histone modifications such as ubiquitination of histone H2B and methylation on histone H3 'Lys-4' (H3K4me3). PAF1C recruits the RNF20/40 E3 ubiquitin-protein ligase complex and the E2 enzyme UBE2A or UBE2B to chromatin which mediate monoubiquitination of 'Lys-120' of histone H2B (H2BK120ub1); UB2A/B-mediated H2B ubiquitination is proposed to be coupled to transcription. PAF1C is involved in mRNA 3' end formation probably through association with cleavage and poly(A) factors. Binds single-stranded DNA (By similarity). Required for maximal induction of heat-shock genes. Required for the trimethylation of histone H3 'Lys-4' (H3K4me3) on genes involved in stem cell pluripotency; this function is synergistic with CXXC1 indicative for an involvement of a SET1 complex.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG215339