

Product datasheet for **MG215376**

Nxf7 (NM_130888) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Nxf7
Synonyms:	nfx- α 1; Nx7; nxf- α ; nxf- α 2; Nxf6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >MG215376 representing NM_130888
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTGCTCTAATGAAAGGAAGCAGGGTCTTCAGATACACAGGAACAATAATGCCTGGCAAAGCCCCA
 TGGGAAGAAAAGAACAGAGGTTTCTCCCTGAGCTATGCTGGCAAGAGGACCGATCATTATGAGCATA
 TGGGCACACGATTCTGCCTTTCAAAGCCCAGAGTGATGGTGGAAATCTGGAGATGAGAGCTGTTCTCTGG
 AACCCGTGAGTAAGACACACCTTATGTTGGGCAACATAGGAGTAGGGAGGAAGAGCATATGGAAGACC
 AAATCCATGTCACGGCGCAGAGAGATGACAAACCTGAACAGAGAGAAGCAGGGCAGTTCACAGAAGACGG
 AGAAGTGGGAGCTGGTTCAAGGTACCGTTCCATGTGGGAGGAAGTATGACAAGACACAGCTAATGAAT
 TCAATCCACAGCCTTTGCAGCGTGCCTTTCACACCTGTGGATTTCCATTGTGACAAAACCGGATCCAGT
 TCTTTGCCCAGACTTTCGATTGCCTCTGCATTGAAGGATATCAGCTATAAGATACATAATGAGTACTT
 TCAGAAGATTCCTATTTTGTCAATCCCTCGTTGCACCTACTCTGTGCAGAATAGATTCACAAAAGAA
 GAAATGGAACAGCTGAAGCTGGCCATGAGGAAAAGATATGATGTTTCCCAGCATGCTCTCTGCCTCAAAA
 AGCTCCGCTTTGATCCAGACTTGATGAAAAACAACATCCACATGATCCTAAATCACAGAAGCTGCATGGC
 TGCCACCTTCAGATCATCCAAGAGAAATTTTCTAAGCTGTTATCCTTGAACCTGAGCAGCAACAACTG
 TTCCAATTAGATAGCCTGATTGATGTGGTAAAGAAGGCACCGCAACTCAAGATCTTGAACCTCCAAAA
 ACATGCTAAGAACAGTGTGGGAGTTGGAGAAGATGAAAGGGCTGAAGCTGGAGCAGCTGTGGCTAGAAGG
 GAACCCATTGTGCAGCACCTTCCCAGACCGATCTTCTATATAAGAGCAGTCTGGAATGTTTTCCGGAG
 TTGTCATACCTGGATGGCAGGAAGTATTACTTCCAAGTGAATGAACACTCAAGAACGCAAGTTAATGA
 AGCCCTGCAAGGACATCTTATGGGTTGAGAGGTGATAAAGAATCAAGTACACAGATTCCTGCGAGAGTA
 TTAAGTGTGATGACTGCTGAAGAGCGACAGGGTCTCCTGAACATTTACCATGACCAGGCTGCTTCTCC
 CTGACCATTCCCTCAACCCCAATGATCCAGATCTGAATAGCATGTATGTCTACTTCAAGGATGAAAATG
 ATATGAAGAATTTCAAGGAGTCCACATTCAGAGGCAGCTACTGAAGTACACAAAACAAGATATTGTGGA
 GTGCTGAGGGGATTTCCCAAACCTACATGCCTTCAAGTCTTCCAAGTGAACATCTGTTTTGAGATG
 GAAACGATGCTCTGCTTCTGTTTGTGGGTTGTTCAAGGAAGAGGGCACCCCAAGGAATGTGTTTCGTG
 CCTTTATGAGGATCTTATTGCTCTCCTTGGCAGTTCAAACATGTATATTGTGAATGACCAGCTTTTTGT
 GAGGAATCCAAGCTCTGAAGAGATCAATAGTGCCTTTGTATACCTCCGCCACATCATACTCCAACCTC
 AAGCTTGCTCTTTCTCAAGAGCAGCAAGAATGGTGCAGGCTTTCTCCACCCAGTCGGGGATGAACTTG
 AATGGTCTCAGAAGTGCCTTGAAGATAACAAGTGGGACTATGCCAGAGCTGCTGAAGTCTTCACTATGCT
 TCAGACCAAGTGAAGATCCCAAAGGAATCTTCAAACAA


ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG215376 representing NM_130888
 Red=Cloning site Green=Tags(s)

MCSNERKQGSFRYTGTNAWQSPMGRKKDRVSSLSYAGKRTDHYEHTGHTILPFKAQSDGGNLEMRVAVLW
 NPSVRHTPYVYGQHRSEEEHMEDQIHVTAQRDDKPEQREAGQFTEDEGELGSWFKVTVPCGRKYDKTQLMN
 SIHSLCSVPFPTVDFHCDKHRIQFFVPDFRIASALKDISYKIHNEYFQKIPFVNPVAPYSVQNRFTKE
 EMEQLKAMRKRYDVSQHALCLKLRFDPDMKNNIHMILNHRSCMAATLQIIQENFPKLLSLNLSNKL
 FQLDSLIDVVKKAPQLKILNL SKNMLRTVWELEKMKGLKLEQLWLEGNPLCSTFPDRSSYIRAVLECFPE
 LSYLDGRKLLLPVMNTQERKLMKPKCKIFMGSEVIKNQVHRFLREYYLMDSEERQGLLNIYHDQACFS
 LTIPFPNDPDLNSMYVYFKDENDMKNFKEFHIQRQLLKYTKQDIVECLRGFPQTLHAFSSFQVNICFQM
 ETMLCFSVCGLFKEEGTPKECVRAFMRIFIALLGSSNMYIVNDQLFVRNPSSEEINSAFVIPSPTSYSNF
 KLVLVSQEQRMVQAFSTQSGMKLEWSQKCLEDNKWDYARAAEVFTMLQTKCKIPKEFFKQ

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:


ACCN: NM_130888

ORF Size: 1860 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_130888.2](#)

RefSeq Size: 1863 bp
 RefSeq ORF: 1863 bp
 Locus ID: 170722
 Cytogenetics: X 57.25 cM

Product images:



Circular map for MG215376