

Product datasheet for **MG215400**

Nolc1 (NM_001039351) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nolc1 (NM_001039351) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Nolc1
Synonyms:	NOPP130; NOPP140; P130
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG215400 representing NM_001039351
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGGATACCGGCTTGCGCCGTGGTTCCAGCGACCTTTATCCCCTTGTGCTCAGATTTCTGCGGG
 ATAGCCAACTCTCGGAGGTGGCCAGTAAATTTGCAAAAAGCGACCGGCGCTACACAGCAGGACGCCAATGC
 CTCGTCCCTCTTGACATCTATAGCTTCTGGCTCAAGTCCACCAAAGCCCCAAAGGTGAAGTTACAGTCA
 AATGGACCAGTGACCAAGAAGGCTAAGAAAAGACTTCATCCAGTGACAGCAGTGAGGACAGCAGTGAGG
 ACGAGGACAAAAAGCCAGGGACTTCCACACAGAAGGCTGCCGCACAGGTCAAGCGAGCCAGTGTGCC
 TCAGCATGCTGGAAGGCGAGCAGCCAAAGCTTCAGAGAGCAGCAGTAGTGAAGAATCCAGTGAGGAAGAG
 GAAGAGGACAAAAAGAAAAGCCTGTCCAGAAGGCAGCTAAGCCCCAAGCCAAGGCAGTCAGACCTCCTG
 CGAAGAAGGCAGAGAGCTCTGAGTCGGACTCAGACTCGGATTCGGACTCCAGCTCAGAGGAAGAAACACC
 ACAGACCCAGAAGCCAAGGCAGCTGTGGCAGCAAAAGCTCAGACTAAAGCCGAAGCCAAACCAGGTACA
 CCAGCGAAAGCACAGCCTAAGGTAGCCAATGGCAAAGCAGCCGCCAGCAGCAGCAGCAGCAGCAGCAGCG
 ATGACTCAGAGGAAGAGAAGAAGGCAGCTGCACCTCCCAAGAAGACTGTACCAAAAAAGCAAGTCGTGGC
 CAAGGCCCCAGTGAAAGTAGCTGCCGCCCCACCAGAAGAGCTCCAGCAGTGAGGATTCTTCCAGTGAA
 GAGGAGGAGGGACAGAGACAACCCATGAAGAAAAAGCAGGTCCCTACAGTTCAGTTCACCACCCTCTG
 TTCCTTTACCAAGAAGTCCCGGGAACCCAGGCTCCAAGAAAGCTGCTGCGCAGACACAGCCTGCAGA
 CAGCAGTGACGACAGCAGTGACGATTCTGATTCAAGTTCAGGAAGAGAAAAACCTCCAGCTAAGACG
 GTCGTCTCCAAGACCCGCCAAAGCAGCTCCAGTGAAGAAGAAAGCAGAAAGCTCTTCAGACAGCTCGG
 ATTCTGACAGTTCTGAGGATGAAGTCTTGCCAAGCCAGTCAGTACAACCAAGAGTCCCAAGCCAGCTGT
 CACTCCGAAGCCATCTGCAGCAAAGGCAGTGACAACCTCTAAGCAACCTGCAGGCAGTAACCAGAAACCT
 CAGAGCAGGAAGGCTGACAGCAGCTCCAGCGAGGAGGAAAGCAGCTCCAGCGAGGAGGAGGAGGCCCTCCA
 AGAAAAGTGCCACAACCCCAAGGCCAAGGTGACTGCTAAAGCAGCACCCGCCAAACAGGCCCTCAGGC
 TGCTGGGGACAGCAGCTCTGACTCAGATAGTTCAGCAGTGAAGAGGAGGAGAAGACTCCTAAGCCCCCA
 GCTAAGAAGAAGGCAGCAGGTGGAGCCGTTTCTACACCAGCCCTGGGAAGAAAGCAGAGGCCAAGAGCA
 GCAGCAGCAGCAGCAGCAGCAGCTCCGAAGATTCCAGTGAAGAGGAGAAAAAAGAAAGCCCAAAGCTAC
 TACCCCTAAAATACAGGCAAGCAAGGCCAATGGCACTCCAGCTTCTCTGAATGGAAAAGCAGCCAAGGAA
 AGTGAGGAGGAAGAGGAGGAGGAAGAAACAGAAGAGAAGAAAAAGGCAGCTGGGACCAAGCCAGGTTTCAG
 GCAAAAAACGGAAGCAGAAATGAGACCGCAGATGAAGCAACAACCTCTCAAGCTAAGAAAGTTAAGCTCGA
 GACCCCAATACGTTTCCAAAAAGGAAGAAGGGAGAAAGAAGGGCGTCTTCCCCTTTCCGAAGGGTCAGG
 GAGGAGGAGATTGAGGTGGACTCTCGAGTGGCGGACAATTCTTTGATGCCAAGCGAGGTGCAGCTGGAG
 ACTGGGGGAGCGAGCCAATCAGGTTCTGAAGTTCACCAAGGCAAGTCTTCCGGCATGAAAAACGAA
 GAAGAAGCGAGGCAGCTACCGGGGAGGCTCCATCTCTGTCCAGGTCAATTCCGTCAAATTCGACAGCGAG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG215400 representing NM_001039351
 Red=Cloning site Green=Tags(s)

MADTGLRRVPSDL YPL VLRFLRDSQLSEVASKFAKATGATQQDANASSLLDIYSFWLKSTKAPKVKLQS
 NGPVTKKAKKETSSSDSSEDSEDEDKKAQGLPTQKAAAQVKRASVPQHAGKAAAKASESSSSEESSEEE
 EEDKKKKPVQKAAKPQAKAVRPPAKKAESESDSDSDSDSSSEETPQTQKPKAAVAQAQTKAEAKPGT
 PAKAQPKVANGKAAASSSSSSDDSEEEKKAAAPPKKTVPKKQVVAKAPVKVAAAAPTQKSSSDSDSSSE
 EEEGQRQPMKKKAGPYSSVPPPSVPLPKKSPGTQAPKAAAQTQPADSSDDSDSDSSEEEKKPPAKT
 VVSKTPAKAAPVKKKAESSSDSDSDSSEDEAPAKPVSTTKSPKPAVTPKPSAAKAVTTPKQPAGSNQKP
 QSRKADSSSEEESSSSEEEEEASKKSATTPKAKVTAKAAPAKQAPQAAGDSSSDSDSSSSEEEETPKPP
 AKKKAAGGAVSTPAPGKKAEEKSSSSSSSSSESDSSEEEKKKPKATTPKIQASKANGTPASLNKAAKE
 SEEEEEEETEKKKAAGTKPGSGKKRQNETADEATTPQAKVKLETPNTFPKRKKGERRASSPFRRVR
 EEEIEVDSRVADNSFDAKRGAGDWGERANQVLKFTKGKSRHEKTKKKRGSYRGGISIVQVNSVKFDSE

TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:


ACCN: NM_001039351

ORF Size: 2100 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039351.2](#), [NP_001034440.1](#)

RefSeq Size: 3646 bp

RefSeq ORF: 2103 bp

Locus ID: 70769

UniProt ID: [E9Q5C9](#)

Cytogenetics: 19 C3

Gene Summary: Nucleolar protein that acts as a regulator of RNA polymerase I by connecting RNA polymerase I with enzymes responsible for ribosomal processing and modification (By similarity). Required for neural crest specification: following monoubiquitination by the BCR(KBTBD8) complex, associates with TCOF1 and acts as a platform to connect RNA polymerase I with enzymes responsible for ribosomal processing and modification, leading to remodel the translational program of differentiating cells in favor of neural crest specification (By similarity). Involved in nucleologenesis, possibly by playing a role in the maintenance of the fundamental structure of the fibrillar center and dense fibrillar component in the nucleolus (PubMed:11424213). It has intrinsic GTPase and ATPase activities (By similarity). [UniProtKB/Swiss-Prot Function]