

Product datasheet for **MG215493**

Slc24a2 (NM_001110240) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc24a2 (NM_001110240) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Slc24a2
Synonyms:	2810021B17Rik; 6330417K15Rik; A1847460; Nckx2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG215493 representing NM_001110240
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATCTCCACCAAAGCCCCACAGCCAGGCTCCTGCAGAAATGGTGTCCCATGAGTCACCATTGGCT
 GCAGGAGACATTATAATAGCAGGAAAAACTGAAGCTAATTCGAGTCATTGGCCTTGTATGGCCCTGGT
 AGCTGTGAGCACTGTCCCATTTTCAATCAGTGCCTTTACTGAGACAGACTCTCAGAGCAACAGAGGAGAA
 GCCAGTGATATGAGTGGCCCCAGAGTAGCACAGGGTCACCGTCAAAGAACTCTCTTAGACTTAAATGACA
 AGATTCGAGATTACCCCCACAGCCACCTGCTTCTCAGGAAGACCAAGCTGAGAACAGTACTGAGCACAC
 CCAGGGGAGACTACCCAAAGGACATCTTTCCCTGGAGGAGCGAAGGAAAGGTGCCATCATTCTGCATGTC
 ATTTGGGATGATCTATATGTTTATAGCATTAGCCATCGTCTGTGATGAGTCTTTGTTCTCTCTAACTG
 TCATCACTGAAAACTGGGAATCTCTGATGATGGCTGGAGCCACCTTCATGGTGCAGGCGGGTCAGC
 CCCAGAACTTTTACGTCTCTTATAGGGGTCTTATTGCTCACAGCAATGTCGGCATTGGCACCATCGTG
 GGTTCGCGGTGTTCAATATTCTCTTTGTTATTGGCATGTGTGCTCTGTTTTCTAGAGAAATCCTAAACC
 TGACCTGGTGGCCGCTCTTCGGGATGTGTCTTCTACATCGTTGATTTGCTCATGTTGATCACATTTTT
 CCTGGATAACGTTATCATGTGGTGGGAAAGCTTACTCCTCTTAACAGCTTATTTGCGCTATGTGGTTTTT
 ATGAAGTTCAACGTCCAAGTAGAGAGATGGGTGAAGCAATGATATCTCGCAATAATGTTATCAAGGTGA
 CAGTGCCAGAAGCCAAAGCAAGCTCTCAGAATCCCTCCAGCCTCACAGGCAGATGGGAGGACAGCTTG
 CCCTACAATAACAAATGACCCGACTGCCAGAGTCTGTGTTTCGAGGTCTACCCTGAGACAGCCAGGCAT
 GGGTCTGCATAACGAATGTGAATGGAGAGCCGCTGAAAGAAAAGACCCAGAACAGCTGTATGGAAGAGT
 CACCTACAGCTGGGGACAAGGATGGGCCAACTTTACCGTGAAGCCACGCCTCAACGAGGCGGAAGCTC
 TGCTCCCTCCACAACAGTCTCATGAGGAATAGCATCTTCCAGCTCATGATACACACCCTTGACCCACTC
 GCTGAAGGGAGGTTAGAGAAAAAGCTTCAATTCTCCACAAGATCGCAAGAAGAAATGCCAAGTGGATG
 AGAATGAGAGGCAGAATGGGGCTGCCAATCATGTGAAAAAATTGAACTCCCAACAGCACAAAGCACAGA
 GGTGGAATGACGCCATCGAGCGAAGCTTCGGAACCCGTCAGAAATGGAACCTCTCCACAACATTGAA
 GCTGCAGATGCCCGAAGGCCACAGAGACTGCCGAGGAAGAAGATGACCAGCCCTCAGCTTGTCTGGC
 CGACCAACACCCGCAAGCAAGCCACCTTCCTAATCGTTTTTCCATTGTGTTCCCTCTCTGGATCACCTT
 ACCTGACGTCGAAAACCTGCTCCAGGAAGTTCTTCCCTATCACATTCTCGGCTCCATCACCTGGATT
 GCAGTTTTCTCATACCTGATGGTGTGGTGGGCACACCAGGTTGGAGAGACAATCGGCATTAGTGAAGAAA
 TCATGGGTCTGACCATCTGGCTGCTGGGACCTCTATCCCTGATCTTATCACCAGTGTATAGTGGCCCG
 AAAGGGCCTGGGGACATGGCTGTGTCCAGCTCTGTTGGAAGCAACATCTTTGACATCACTGTGGGGCTC
 CCACTGCCCTGGCTCCTACACCATCATCCACAGATTAGCCAGTACTGTCAGCAGCAATGGCTCT
 TCTGTGCCATCGTCTTCTTTCATCATGTTGCTTTTGTCTTCTCCATCGCCCTCTGCAAATGGCG
 GATGAACAAAATCCTCGGCTTATCATGTTTGGCCTCTACTTCGTGTTCTGGTGGTCCAGGTCCTCCTG
 GAAGACAAAGTTCTCGTGTGCCCTGTCTCCATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG215493 representing NM_001110240
 Red=Cloning site Green=Tags(s)

```
MDLHQSPARLLQKWCSESHPFGCRRHYNSRKKLKLIRVIGLVMGLVAVSTVPFISAFETETDSQSNRGE
ASDMSGPRVAQGHQRQRTLLDLNDKIRDYTPQPASQEDQAENSTEHTQGDYPKDIFSLEERRKGAILHV
IGMIYMFIALAIVCDEFFVPSLTVITEKLGISDDVAGATFMAAGGSAPELFTSLIGVFIHNSVIGITIV
GSAVFNILFVIGMCALFSREILNLTWWPLFRDVSFYIVDLLMLITFFLDNVIMWESLLLLTAYFAYVVF
MKFNVQVERWVKQMISRNNVIKVTYPEAQAKLSEFPASQADGRTACPTITNDPDCQSLCFEVYPETARH
GSCITNVNGEPLKEKTQNSCMEESPTAGDKDGPTLPSKPRLQRGSSASLHNSLMRNSIFQLMIHTLDPL
AEGRFREKASILHKIAKKKQVDENERQNGAANHVEKIELPNSTSTEVENTPSSEASEPVQNGLSHNIE
AADAPKATETAEEDDQPLSLSWPTNTRKQATFLIVFPIVPLWITLDPVRKPARKFFPITFFGSIWI
AVFSYLMVWWAHQVGETIGISEEIMGLTILAAGTSIPDLITSVIVARKGLGDMAVSSVSGNSIFDITVGL
PLPWLLYTIHRFSPVTVSSNGLFCAIVLLFIMLLFVILSIALCKWRMKNILGFIMFGLYFVFLVSVLL
EDKVLVCPVSI
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001110240

ORF Size: 2133 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001110240.1, NP_001103710.1


RefSeq Size: 10516 bp

RefSeq ORF: 2136 bp

Locus ID: 76376

Cytogenetics: 4 C4

Product images:



Circular map for MG215493