

Product datasheet for **MG215494**

Uqcc3 (NM_001160356) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: TurboGFP

Symbol: Uqcc3

Synonyms: AI462493; Cbp4

Mammalian Cell Selection: Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG215494 representing NM_001160356
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGGTGGCTCGTAAAGCACTTGTGGCAGTTGCACTGCTAGGCGGGGAGCTGGCGTGGTTCTATTC
TGTTTGTCTTGTGACCCAGGAGAACTACAGAAGCAGTCGATGCTGCAGGAGATGCCGAAAGGGACTC
GCGGCGCAGGGACGAAGCAGTCAGGACCACGGAAGTGGTGTGGCTACCCTGAAGGACGCCGACGCCACG
AAGGAGAACGTGGCCTGGAGGAGAACTGGACAGTTAGCGGGGATGGCAGGTCAGCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG215494 representing NM_001160356
Red=Cloning site Green=Tags(s)

MEVARKALVAVAVLGGGAGVGSILFALVTPGELQKQSMQEMPERDSRRRDEAVRTELVMATLKDAAAT
KENVAWRRNWTVSGDGRSA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: NM_001160356

ORF Size: 267 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001160356.1 , NP_001153828.1
RefSeq Size:	847 bp
RefSeq ORF:	270 bp
Locus ID:	107197
UniProt ID:	Q8K2T4
Cytogenetics:	19 A

Gene Summary: Required for the assembly of the ubiquinol-cytochrome c reductase complex (mitochondrial respiratory chain complex III or cytochrome b-c1 complex), mediating cytochrome b recruitment and probably stabilization within the complex. Thereby, plays an important role in ATP production by mitochondria. Cardiolipin-binding protein, it may also control the cardiolipin composition of mitochondria membranes and their morphology.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG215494