

Product datasheet for **MG215516**

Olfr376 (NM_001172686) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr376 (NM_001172686) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Olfr376
Synonyms:	MOR135-12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG215516 representing NM_001172686 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAGCCAGATGCAGGGCACCACAGAAAGGAATCAAAGTCCATCTCCAGTTCCTCCTCTGGCCTGC
CCATCCCCACAGAGCATCAACACCTGTTCTATGCCCTGTTCTGGCCATGTACCTCACCACCGTCTGGG
GAACCTCATCATCATCCTTATTCATCTGGACTCCCATCTCCACACCCATGTATTCTTCTCAGC
AATTTGCTCTCTGATCTCTGCTTTTCATCTGTGCAATGCCAAATTGCTGCAGAACATGCAGAGCC
AGGATCCATCAATTCCCTATGCAGGCTGCCTGGCACAATGTACTTTTCTGTTTTTGCAGACCTTGA
AAGCTTCTCCTTGTGGCCATGGCCTATGACCGCTATGTGGCTATCTGCTTCCCCCTTCATTACATGAGC
ATCATGAGCCCCAGGCTCTGTGTGAGTCTGGTGTCTGTCTGGGTGCTGACCACCTTCCATGCCATGC
TGCATACCCTGCTCATGGCCAGATTGTCTTGTGTGAGGACAATGTATCCCCACTTTTTCTGTGACAT
GTCTGCTCTGTGAAGCTGTCTGCTCTGACACCTATGTTAATGAATTGGTGATATTTGTGATGGGAAGC
CTGATCCTTGTGATTCCATTTGTGCTCATCCTTGTGCTATGCACGAATTGTGCTCTCCATTCTCAAGG
TCCCGTCTGCTCGAGGCATCCGTAAGCCTTCTCCACCTGTGGGTCCCACCTGTCTGTGGTGTCACTGTT
CTATGGAACAGTCATTGGTCTGTACTTATGTCCATCAGCTGATAACTCTACTGTGAAGGAACTGTCATG
GCCATGATGTACACAGTGGTGACTCCCATGCTGAACCCCTTCATCTACAGCCTGAGGAACAGAGACATGA
AGGGGGCCCTAATTAGTGTCTTTGCAAGAAGAAATTCTTTTCTGTCTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG215516 representing NM_001172686
 Red=Cloning site Green=Tags(s)

MSQMGGTTERNQTAISQFLLLGLPIPTHEQHLYALFLAMYLTTVLGNLIIIIHLDSHLHTPMYSFLS
 NLSFSDLCFSSVTMPKLLQNMQSQDPSIPYAGCLAQMYFFLFFADLESFLLVAMAYDRYVAICFPLHYMS
 IMSPRLCVSLVLLSWLTTFHAMLHTLLMARLSFCEDNVIPHFCDMSALLKLSCSDTYVNELVIFVMGS
 LILVIPFVLILVSYARIVSSILKVPVSARGIRKAFSTCGSHLSVSVSLFYGTVIGLYLCPADNSTVKETVM
 AMMYTVVTPMLNPFIIYSLRNRMKGALISVLCCKKILFCL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001172686

ORF Size: 960 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001172686.1, NP_001166157.1

RefSeq Size: 1210 bp


RefSeq ORF: 963 bp

Locus ID: 258924

Cytogenetics: 11 B4

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG215516