

## Product datasheet for **MG215519**

### **Olf52 (NM\_146583) Mouse Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Olf52
Synonyms:	IE6; MOR185-6; ORL446
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG215519 representing NM\_146583  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCTCAGATAAACTGCACCCAGGTGACAGAGTTTATTCTTGTTGGCCTCACAGATCGCAAGGAGTTAA  
AGATGCCCCCTTTGTGGTGTCTTATTTCATCTATCTTTTCACAGCAATAGGCAATCTGGGATTGATCTT  
GGTCATTAGAACAGATGCAAGACTCAACACACCCATGACTTCTTCTCAGCAACCTGGCTTTTGTGAT  
TTCTGCTACTCCTCTGTCACTACACCTAAGATGCTTGGGAATTTCTGTATAGCAAAAATGCCATATCCT  
TCAATGCATGTGCTGCCAGTTAGGCTGCTTTCTCACATTTATGGTATCAGAGTGCTTGTCTCTGGCTTC  
CATGGCATATGATAGATATGCAGCAATTTGTAACCCCTCTATTGTATATGGTCACAATGTCTCCTGGAATC  
TGCATTAGCTTGTAGTTGTGCCCTATAGCTATAGTTTCTCATGGCATTGATTCACACTCTTCTAACCT  
TCCGCCTATCCTATTGCCATTCTAATATCATCAATCACTTCTACTGTGATGACATGCCTCTTCTCAGGCT  
AACTTGCTCAGATACTCACTACAAGCAGCTGTCTATTTTGGCCTGTGCTGGAATCACATTCTTTCTTCT  
GTTCTGATTGTTTCTGTATCCTACATGTTTCATTATTTCTGCCATTCTGAGGATGCGCTCAGCTGAAGGAA  
GACGAAAGCCTTTTCCACCTGTAGCTCTCACATGATGGCAGTGAGCATATTCTATGGAACCTTATCTT  
TATGTACTTACAGCCGAGCTCTGACCATTCTTTGATACAGATAAGATGGCCTCTGTCTTCTACACAGTG  
ATCATCCCATGTTGAACCCCTTGTCTACAGCCTCAGGAACAAGGATGTGAAAGATGCCCTGAAGAGGG  
TCATGGACAATAGAAACAGACTTCTATATTTAGAGAACTAAGAAAA

**ACCGGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



**Protein Sequence:** >MG215519 representing NM\_146583  
Red=Cloning site Green=Tags(s)

MAQINCTQVTEFILVGLTDRKELKMPLFVVFLFIYLFITAIGNLGLILVIRTDARLNTMPYFFLSNLA FVD  
 FCYSSVITPKMLGNFLYSKNAISFNACAAQLGCF LTFMVSECLLLASMAYDRYAAICNPLLYMVTMSPGI  
 CIQLVVPYSYFLMALIHTLLTFRLSYCHSNIINHFY CDDMPLLRLTCS DTHYKQLSILACAGITFISS  
 VLIVSVSYMFIISAILMRSAEGRKAFSTCSSHMMAVSIF YGTLIFMYLQPS SDHSLD TDKMASVFYTV  
 IIPMLNPLIYSLRNKDVKDALKRVM DNRNQTSIFREL RK

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**ACCN:** NM\_146583

**ORF Size:** 957 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_146583.1](#), [NP\\_666794.1](#)

**RefSeq Size:** 960 bp


**RefSeq ORF:** 960 bp

**Locus ID:** 18352

**Cytogenetics:** 2 D

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG215519