

Product datasheet for **MG215525**

Olf1384 (NM_146472) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Olf1384
Synonyms:	MOR256-23; RP23-58E13.3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG215525 representing NM_146472
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**


ATGGGAAGTTTCAATCTCAGCTTAGAAGGTTTCTGTTGGTAGGCTTCTCAGACTGGCCTCAACTGGAAC
TTATCCTCCTGGTCTTCATTTTCGATTTTCTACTCTTACCCTCTGTGGCAACATCACCATTATTGTTCT
CACTCAACAAGACCTTCATCTGCACACACCCATGTACTTCTTCCCTGCCCACCTCTCCTTCTGGACCTC
TGCTTCACCAGTAGCACTGTGCCAAACTTCTAATCAGTCTTCCCGAGGTGATCAGACTATCAGCTATG
CTGGGTGCATGACCCAGTTCATAGCACTTCTACTGGGTGGAAGTGTGTGTGCTCCTTGTGGTGAT
GGCTTTCGACCGCTATGTCGCTGTGTGTCGCTCCACTACACTACACCAGCATTATGCACCCCTTCTCTGC
CATGCATTGGCCATCTCCTCCTGGGTGGGAGGCCTTGTGAAGTCTCTGACTCAGACAAGTCTCATCATGA
CCATACCTCTCTGTGGCCATCACCTGAACCACTTCTTCTGTGAGATGCTTGTCTCCTGAAGCTGGCTTG
TGAGGACACAGTGGGAACAGAGACCTACTTGTGTTGGGCTGGAGCTGTAATTTTGGTCTGTCTGTAGCA
CTAATTCTAGGCACCTATGCACACATTGCTCATGCAGTGTGAAGTCAAGTCAAGTCTGGGCGCAGAA
AGGCTCTGGAACTTGTGGTCCCACCTTACTGTGTTTTTCTTTTTATGGCTCAGCCATGTACACGTA
TCTCCAGCCTATCCACACGTATTCTGGTAGTGAAGGAAGTTTGTGCCCTCTTTTATACTATCATTACT
CCAATGCTGAACCTCTGATTTATACCCTAAGAAACAAGGATGTAAGGGGGCTCTGTGCAAGGACTTGT
TAAAGGAAAAAAGGAGACTAAGACTAGAAGAAAGATGTTGAA

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_146472.2</u> , <u>NP_666683.1</u>
RefSeq Size:	1034 bp
RefSeq ORF:	957 bp
Locus ID:	258464
Cytogenetics:	11 B1.2
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG215525