

Product datasheet for **MG215572**

Olfr159 (NM_019476) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr159 (NM_019476) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Olfr159
Synonyms:	mOR37e; MOR262-5; Olfr37e
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG215572 representing NM_019476 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGAGATCCAACAAGACCACCCCTGTGTCCAGCTTCATCCTCCTGGGTCTCTCCGCCACCCAAAGC
 TGGAGAAAACCTTCTTCGTGCTCATCCTGCTGATGTACCTGGTGATCCTGCTGGGCAATGTCGTCTCAT
 CCTGGTGAGCATCCTCGACTCCACCTGCACACGCCCATGTACTTCTCCTGGGGAACCTCTCCTTCCTG
 GACATCTGCTACACCACCTCCTCCGTCCCCCTCATTCTGGACAGCTTCTGACCCCCAGGAAGACCATCT
 CCTTCTCGGGCTGTGCCGTGCAGATGTTTCTCTCCTTCGCCATGGGAGCCACAGAGTGTGTCTCCTGGG
 CATGATGGCGTTTGATCGTTATGTGGCCATCTGCAACCCCTTAGATATCCTGTGGTCATGAGCAAGGCT
 GCCTATGTGCCCATGGCCGCCGGCTCCTGGGTATCTGGAAGCATCACTGCCACAGTGCAGATATCTTTAG
 CAATGACACTGCCCTTCTGTGGGACAATGTCATCAATCACTTTACATGTGAGATCCTGGCAGTCTGAA
 ACTGGCCTGTGCTGACATCTCCATCAATGTCATCAGCATGGCAGTGGCTAATGCAATGTTTTGGGGGTC
 CCTGTCCTTTCATCTTTGTCTCCTACATCTTCATACTCTCCACCATCCTGAGGATCCCCCTGCGGAGG
 GGAGGAAGAAGGCCTTCTCCACCTGCTCTGCCACCTACCGTGGTACTTGTCTTCTATGGGACCATCCT
 GTTCATGTACGGGAAGCCAAATCCAAGGACCCATTGGGGGCAGACAAACAGGACCTTGCAGACAAGCTC
 ATCTCCCTCTTCTATGGAGTGGTGACCCCATGCTGAACCCATCATCTACAGCCTGAGGAACAAGGATG
 TGAAGGCCGCTGTGACAAACCTGGTGGCCAGAAACACTTCAAGTGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >MG215572 representing NM_019476
 Red=Cloning site Green=Tags(s)

MERSNKTPVSSFILLGLSAHPKLEKTFVLLIMYLIVLLGNVLLILVSILDSHLHTPMYFFLGNSFL
 DICYTSSVPLILDSFLTPRKTSISFGCAVQMFLSFAMGATECVLLGMMAFDRYVAICNPLRYPVMSKA
 AYVPMAGSWVSGSITATVQISLMTLPFCGDNVINHFTCEILAVLKLACADISINVISMAVANAMFLGV
 PVLFI FVS YIFILSTILRIPSAEGRKKAFSTCSAHLTVVLVFGYGTILFMYGKPKSKDPLGADKQDLADKL
 ISLFYGVVTPMLNPPIIYSLRNKDVKAAVTNLVGQKHFKW

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_019476

ORF Size: 957 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_019476.1, NP_062349.1

RefSeq Size: 960 bp


RefSeq ORF: 960 bp

Locus ID: 29849

Cytogenetics: 4 A5

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG215572