

Product datasheet for **MG215641**

Crip3 (NM_053250) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Crip3 (NM_053250) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Crip3
Synonyms: TImp; TLP; TLP-A; TLP-B
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG215641 representing NM_053250
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGCTGGACTTGTCCGCGTTGTCAGCAACCCGTTTACTTTGCTGAGAAGGTGAGCTCCTTAGCAAGA
 ACTGGCACCCTTCTGCCTGAAATGTGAGCGCTGTACAGCATCTTGTCCCCTGGCGGCATGCAGAGCA
 CAACGGGAGACCATACTGCCACAAGCCTTGCTATGGGGCTCTTTGGACCTAGAGGGGTCAACATTGGT
 GGTGTGGGCTGCTATCTCTACAACCTCCCAACTCCCCCTCTGCCAGCCGATTTCCCTCAGCCCCAGCA
 ACTTTAGCCCCCAGGCAAGGACTGGCCTCTCCCGGCCAAGAAAAGGCTCCCTACTGAAGACATT
 CACTGGTGAGACTTCACTGTGCCCTGGATGTGGGGATCCTGTCTTTTTGCTGAGAAAGTGATGTCTTTG
 GGCAGAAATTGGCACCGACCCTGCCTGAGGTGCCAACGTTGCAGGAAGACCCTGACTGCTGGGAGTCATG
 CGGAGCATGATGGCATGCCTTACTGCCACATCCCCTGCTATGGATACCTGTTTGGCCCCAAGGTGGGCA
 ATCTGACCCAGACACTGGACCTATGAGAATGTGTGGACAACCTCCCTTGTGGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG215641 representing NM_053250
 Red=Cloning site Green=Tags(s)

MSWTCPRCQQPVYFAEKVSSLGKNWHRFCLKCERCHSILSPGGHAEHNGRPYCHKPCYGALFGPRGVNIG
 GVGCYLYNLPPTPPASRISLSPSNFSPPRPRTGLSRAKKRPPYLKFTTGETSLCPGCGDPVFFAEKVMSL
 GRNWHRPCLRCQRCRKTLAGSHAHDGMPYCHIPCYGYLFGPKGGQSDPRHWTYENVWTTSLVV

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_053250

ORF Size: 615 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053250.2](#), [NP_444480.1](#)

RefSeq Size: 1014 bp


RefSeq ORF: 618 bp

Locus ID: 114570

UniProt ID: [Q6Q6R3](#)

Cytogenetics: 17 C

Product images:



Circular map for MG215641