

Product datasheet for **MG215700**

Mak16 (NM_026453) Mouse Tagged ORF Clone

Product data:

| | |
|--------------------|--------------------------------|
| Product Type: | Expression Plasmids |
| Tag: | TurboGFP |
| Symbol: | Mak16 |
| Synonyms: | 2600016B03Rik; AI31491I; Rbm13 |
| Mammalian Cell | Neomycin |
| Selection: | |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |

ORF Nucleotide Sequence: >MG215700 representing NM_026453
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCAGTCGGACGACGTAATCTGGAATACACTAGGAAACAAGCAGTTTTGTTCTTTCAAATAAGAACCA
AGACTCAGGGCTTCTGCAGAAATGAATACAGTCTGACTGGTCTGTGCAATCGGTCTGCTCTGCTCTGGC
AAATAGTCAGTATGCTACCATTAAGAAGAGAAAGGGCAGTGTACCTGTACATGAAGGTTATAGAGCGA
GCAGCTTTTCCAGGCGTCTCTGGGAACGGTCCGGCTTAGTAAAACTATGAAAAAGCACTGGAACAAA
TAGATGAAAACTGATCTATTGGCCCGCTTATCAGACACAAATGTAAGCAGAGGTTTACCAAGATAAC
CCAGTACCTGATTGGAATCAGAAACTTACACTGAAGCGACAGAGGAAGCTTGTTCCTTTGAGTAAGAAG
GTTGAACGAAGGAAAAACGAAGAGAGGAAAAAGCATTAAATAGCTGCCAGCTGGACAATGCTATTGAGA
AGGAATTGCTGGAGAGACTGAAACAAGATACGTATGGCGACATCTACAACCTCCCCATCCATGCCTTCGA
CAAGGCCCTAGAGAAACAGGAAGCAGAAAGTGACTCTGAAAGATGAAGAAGAAGAAGATGAAGACGAA
GAAGAGGATGTGGGAAAAGAGAGTTTGTAGAAGATGAGGAGTGGAGGAGAGTGACCTGAGTGACTTTG
AGGATATGGATAAACTGAATACTGACAGTGAGGAAGACCAGGATGATGAATCCTCCAATGATGAAGAGGC
TCATAAGGCTAAGCACAAAGGCAAAGCACCCCTTGAAAGGACCTTACGGAAAAACGGGCCATATGTGAA
ATAGAGTATGAACAGGACAGAGCCCATGGCCAAAGTCAAAGCCACA

ACCGGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >MG215700 representing NM_026453
Red=Cloning site Green=Tags(s)

MQSDDDIWNTLGNKQFCFSFKIRTKTQGFRCNEYSLTGLCNRSSCLANSQYATIKEEKGCQCYLYMKVIER
 AAFPRRLWERVRLSKNYEKALEQIDENLIYWPRFIRHKCKQRFTKITQYLIRIRKLTLLKRQRKLVPLSKK
 VERREKRREKALIAAQLDNAIEKELLERLKQDITYGDIYNFPIHAFDKALEKQEAESDSEDEEEEEDEDE
 EEDVGKREFVEDEEVEESDLSDFEDMDKLNNTDSEEDQDESSNDEEAHKAKHKGKAPLKGPLRKKRAYVE
 IEYEQETEPMAKVKAT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:




ACCN: NM_026453

ORF Size: 888 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| Note: | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required. |
| RefSeq: | NM_026453.3 , NP_080729.1 |
| RefSeq Size: | 1857 bp |
| RefSeq ORF: | 891 bp |
| Locus ID: | 67920 |
| UniProt ID: | Q8BGS0 |
| Cytogenetics: | 8 A3 |

Product images:

Circular map for MG215700