

Product datasheet for **MG215761**

Taar8b (NM_001010837) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Taar8b
Synonyms:	Gm1149; taR-8b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG215761 representing NM_001010837
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACCAGCAACTTTCCCAACCAGCCCTGCAGCTCTGCTATGAGAACACGAATGGATCCTGTATTA
CTCCCTATTCTCCAGGGCCCCGGGTCACTCTGTACATGGTCTTTGGCTTTGGGGCTGTGCTGGCAGTGTG
TGGGAACCTCCTGGTGGTGATTTTCAGTTCTCCATTTCAAGCAGCTGCACTCTCCAGCCAATTTCTCATC
GCCTCTCTGGCCAGCGCTGACTTCTTGGTGGGCATCTCTGTGATGCCCTTCAGCATGGTCAGGTCCATTG
AGAGCTGCTGGTACTTTGGAGATGCATTTTGTAGCCTTCACAGTTGCTGTGATGTGGCTTTTTGTTACTC
TTCTGCCCTCCACCTGTGCTTCATCTCAGTGGACAGGTACATTGCTGTCACAGACCCTCTGGTCTATCCC
ACCAAGTTCACAGTGTCTGTCTGGAATTTGCATCAGCATCTCCTGGATTCTGCCCTGGTATACAGCA
GTGCAGTGTCTACACAGGCATCAGTGTCTAAGGGATTGAAAGCTTAGTAAGTCTCTGAATTGTGTAGG
GGGCTGCCAAATCGTTGTCAATCAAGACTGGGTTTTGATAGATTTTCTTCTATTTTTCATACCTACCCCT
GTTATGATCATTCTTTACAGCAAAATTTTTGGTAGCCAAACAGCAAGCTGAAAAATTGAAACTTCTG
TAAGTGACAACAGAGGTGAATCATCTCAGAGAGTCACAAAGCCAGAGTGGCCAAGAGAGAGAGGAAGGC
TGCAAAAACCTTGGGGTCACTGTGGTGGCTTTTATGGTCTCGTGGCTCCCATACACAATTGACTCATTG
GTTGATGCTTTTCGTGGCTTCATTACTCCTGCCTATGCTATGAAATCTGTTGCTGGAGTGCCTATTATA
ACTCAGCCATGAACCCCTTTGATTTATGCTTTCTTTTATCCTTGGTTTTAGGAAAGCCATAAAGCTTATTTT
AAGTGGGGAGATTCTAAAGAGTCATTCATCTACCATGAGTTTGTTCAGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >MG215761 representing NM_001010837
Red=Cloning site Green=Tags(s)

MTSNFSQPALQLCYENTNGSCIKTPYSPGPRVILYMFVFGAVLAVCGNLLVVISVLHFKQLHSPANFLI
 ASLASADFLVGISVMPFSMVRSDIESCWYFGDAFCSLHSCCDVAFYSSALHLCFISVDRYIAVTDPLVYP
 TKFTVSVSGICISISWILPLVYSSAVFYTGISAKGIESLVSALNCVGGCQIVVNQDWLIDFLFFIPTL
 VMIIILYSKIFLVAQQAVKIETSVSDNRGESSSESHKARVAKRERKAAKTLGVTVVAFMVSWLPYTI
 VDAFVGFITPAYVVEICCSAYNSAMNPLIYAFFYPWFRKAIKILSGEILKSHSSTMSLFSE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001010837

ORF Size: 1032 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001010837.1 , NP_001010837.1
RefSeq Size:	1035 bp
RefSeq ORF:	1035 bp
Locus ID:	382348
UniProt ID:	Q5QD06
Cytogenetics:	10 A4
Gene Summary:	Orphan olfactory receptor specific for trace amines.[UniProtKB/Swiss-Prot Function]

Product images:

Circular map for MG215761