

## Product datasheet for **MG215812**

### Gsx1 (NM\_008178) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gsx1 (NM_008178) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gsx1
Synonyms:	Gsh-1; Gsh1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG215812 representing NM_008178 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCGCGCTCCTTCTGGTGGATTCCCTTGTGCTGCGGGAAGCCAGCGACAAGAAGGCTCCGGAGGGCA  
GCCCGCCACCGCTCTTCCCCTACGCGGTCCCGCCCGCCAGCGCTCCACGGCTCTCGCCGGGCGCCTG  
CCACGCGCGCAAGGCCGGCTTGTGTGCGTGTGTCCCCTCTGTGTACCAGCTTCGCAGCTGCACGGGCC  
CCCGGGCCGCGGCACTGCCGCTACTCAAGGCGTCCCTCCCTCCGATCGCAGTACTGCCACGCAC  
CCCTGGGCGCCAGCACTCCGTGTCCCCTGGAGTCGCCACGGCCCGGCCGCGGCCGAGCAGCAGCTGC  
ACTCTACCAGACCTCCTACCCGCTGCCGATCCCAGACAGTTTCACTGCATCTCTGTGGACAGCAGCTCG  
AACCAGCTGCCAGCAGCAAGAGGATGCGGACGGGTTACCAGCACACAGCTCCTGGAGCTGGAGCGAG  
AGTTTCGCCTCCAACATGTACCTCTCCCGCTGCGGCGCATCGAGATCGCGACCTATCTGAACCTGTCCGA  
GAAGCAGGTGAAGATCTGGTTTCAGAACC CGCGGTGAAGCACAAGAAAGGCAAAGGCAGTAACAC  
CGCGGCGGAGCTGGGGCGGGGCGCGGGGGCGCACCGCAAGGCTGCAAGTGCTTTCGCTCTCCTCAG  
CCAAATGCTCAGAGGACGACGACGAATTGCCATGTCTCCATCTTCTCCGGGAAGGATGACAGAGATCT  
CACAGTCACTCCG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG215812 representing NM\_008178  
 Red=Cloning site Green=Tags(s)


MPRSFLVDSLVLREASDKKAPEGSPPPLFPYAVPPPHALHGLSPGACHARKAGLLCVCPLCVTASQLHGP  
 PGPPALPLLKASFPPFGSQYCHAPLGRQHSVSPGVAHGPAAAAAAAAAALYQTSYPLPDRQFHCISVDSSS  
 NQLPSSKRMRTAFTSTQLLELEREFASNMYLSRLRRIEIATYLNLSKQVKIWFQNRVRVHKHKEGKGSNH  
 RGGAGAGAGGGAPQGCKCSSLSSAKCEDDELPMSPSSSGKDDRDLTVTP

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_008178

**ORF Size:** 783 bp

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_008178.2</a> , <a href="#">NP_032204.1</a>
<b>RefSeq Size:</b>	1773 bp
<b>RefSeq ORF:</b>	786 bp
<b>Locus ID:</b>	14842
<b>UniProt ID:</b>	<a href="#">P31315</a>
<b>Cytogenetics:</b>	5 86.79 cM
<b>Gene Summary:</b>	Probable transcription factor that binds to the DNA sequence 5'-GC[TA][AC]ATTA[GA]-3'. Activates the transcription of the GHRH gene. Plays an important role in pituitary development.[UniProtKB/Swiss-Prot Function]