

Product datasheet for **MG215847**

Lrrc36 (NM_001033371) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc36 (NM_001033371) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Lrrc36
Synonyms:	BC036564
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG215847 representing NM_001033371
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCTCACGCCGTGGCTTTCGCCTCCAGATGACCGAGCCGTCCGTGACAGCGAGAGGAGAGCGGCCAAGC
 TACATTTTCAGTCAGTTGGGCAACAGTGAATTTTCTCTTAGAGGTGAAAAAAGCTCTAGGGAGAAGAC
 AATGAAAACTGTGTGACAGATGAGGGCCCGCTCACACGTCACTCCTGAGTTGATGCCAGGATGGAA
 ACAGACGCAAAACAAAGGACTTTTTATTCCTTCCCAACCGGAAATAAAGGATTCCCTAACAGCATTT
 GTGCAGCCAGGGCAGTGGCACACCAGCTCAGAAATTAGACGTTTTTCCGCTGGGACACAGATGCAGGA
 AGCAACAAGAAGGGAGACAAGTGACCTTACCAGGAAGATGGATTAGAACTCAGGCTCTACCCGCTCCT
 CAGAGCACAGTCCGATCTCCAGAGAAGATGACGAGAGATGGATACCGAGTGTCTTCTAGACACTAAGT
 CTTCAGGTTCTTCTCCAGAAAAGGACTTGATACCAAAACCTGATGCCTATCCATGTACCCATGATGCCTC
 CTTGGGTAAACGCTGGATGTGGGTGATTCTAACCAGATCCTTCCCTGTCAGTTACCTTCTGAGGTTTGT
 CTGGACCACTATGGTAATCAATATTCTCAAACCTTTGCCTGCACGGAAGCTTGTCAAAAGGGCTCAGA
 AGGGCAAGAACTACCGAGAACATAGCATAAAGCCTTACAGGACAAAAAGGCCACCATCCCATCCCTG
 TGGAGACTTAAACATCCCTGTCAAACCCGACTCCAGCACGGGAAGGCTGCTGAGGCTTAGCTCAGAT
 CTCTATGCCACAACCCATTTCAACAGTGACCCCGCTCTGCTGGCAACGTGGAGCAACAGCTTTCAGCC
 TCCGTGACTTTACACCAGCACCTGGTCTTTCCCAAGCAGCCCTGCCCTAGGAAATTCTCTGCGGACACT
 GCTGCTGCCTCCTGGGACTCCAGAAAACAGGGAGATTCGACTAAGAGGTCATTAAGCCCTTCAAGGAGA
 GGATTCAAGCGGAAGGACGGTATCCTTGCAACCCGAGTCTGAAGCGTGGTTTCAAGATGCTACAGGCA
 GCGAGGCGCAACCTCTCTAGTGACCTGGGCACTTTGCATGGTTATCAGGAAACACAGTCTCCGAT
 CTCTGCCAGAACCCCGCTGGCCACTGTCCTCAGACAGCTCCTGGAGCTGGTGGACAAGCACTGGAAT
 GGCTCTGGCTCCCTCCTCTGGACAAGAAGTTTCTCGGTCCCGCCAGAGACTTGTCTCTGTCATTGGTGG
 TCCCTGCACCCTCTCAGCAGTGGCGGCGCTCGAAGCTTGACGATAAGGCAGGGAAGCCTTGTGCTGGAG
 GGAGACCGAGCTAAAGGAAGCAGGGCTGCTAGTTCCTAACGATGTGGAAGTCTGAAGCAAAACTGGTC
 AAAGTGTGGAGGAAAACCTTGTCTGTGCAAAAAATTCAGCAGTTGGAGGGAACCGCTGCCACCTCGA
 TTGTGAGTGGGACCCATCACACACACGATGAGCTTCTGCGCAAAACCAACAGCTGACCATACAGGT
 GGCTTGCCTGACCCAGGAGCTGACTCAGCTGAAAAGACTGGAGGAGACAGTTGCTCTTCTCCATGAGAGT
 CAGAGGTCCCTGGTGGTAACTAATGAGTATCTGCTGCAGCAGCTGCATAAAGAGCAGAAAGTTATCCG
 GGAAATCTCTCCTGCCTCCCGAGAAGAGTCATCTTTGGGGAGATCGTCGCCCTTTGGGAAAAGCAGCT
 GTCTTCTCCTCACCATGGTACACGACACGGGTCACTACCTAATACAGAGTGTCTCAGAAGCCGACCCA
 GAGCTAGTTTGTGGAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:


>MG215847 representing NM_001033371
 Red=Cloning site Green=Tags(s)

MLTPWLSPPDDRAVRDSERRAAKLHFSQLGNSENFLLVEKSSREKTMKNCVTDEGPASHVSPEVDARME
 TDANKGLFIPFPNREIKDSLTSICAAQSGTTPAQKLDVFPLGTQMGEATRRRSDLHQEDGLELRLYPPP
 QSTVRSPEKMTDGYRVSFLDTKSSGSSPEKDLIPKPDAYPCTHDASLGKRLDVGDSNQILPCQLPSEVC
 LDHYGNQYSQTLCLHGSLVKRAQKGKGYREHSIKPSQDKKATTSHPCGDLLTSLNPDSSSTRLLRLSSD
 LYATTHFNSDPALLANVEQQLSLRDFTAPGSPFPSSPALGNSLRTLLLPPGTENREIPTKRSLSPSRR
 GFKRKDILANPSLKRGFQDATGSEAQPLSSDLGSLHGLSGNHSPPI SARTPHVATVLRQLLELVKHWN
 GSGSLLLDKKFLGPARDLLLSLVVPAPSQQWRRSKLDDKAGKALCWRETELKEAGLLVPNDVESLKQKL
 KVLEENLVLEKIQQLLEGTAATSIIVSGHPSHTHDELLRKNQQLTIQVACLTQELTQLKRLEETVALLHES
 QRSLVVTNEYLLQLHKEQKGYSGKSLLPPEKSHPLGRSSPFKSTLSSSSPMVHDTGQYLIQSVSEADP
 EPSLWS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001033371

ORF Size: 1908 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001033371.3](#), [NP_001028543.1](#)

RefSeq Size: 2102 bp


RefSeq ORF: 1905 bp

Locus ID: 270091

UniProt ID: [Q3V0M2](#)

Cytogenetics: 8 D3

Product images:



Circular map for MG215847