

## Product datasheet for **MG216124**

### Tmem158 (NM\_001002267) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tmem158 (NM\_001002267) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Tmem158  
**Synonyms:** 2310037P21Rik; Ris1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG216124 representing NM\_001002267  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGCTGCCCTTCTCGCCGCGCTACTGGCCGCGCCTGCCAGCTGCCGCCGGCGCACGGCGGGGCCACGG  
 ATGCGCCCGCCTTGCAGGGACGCCCCAAACGCCTCGGCCAATGCGTCCTTCACCAACGAACACTCCAC  
 CCCGCTCTGCTGGCATCGGCGGCATCGGCGCCCCAGAGCGCGCGGGCCAGAGGAGGCTCCGGCGGGC  
 CCCTGCAACATCAGCGTGCAACGACAGATGCTGAGCTCGCTGCTCGTGCCTGGGGCCGCCCGGGGCT  
 TGCAGTGCACCTGCTGCTCTTCTCTACCAACGCGCACGGCCGCGCCTTCTTCGCCGCGCCTTCCACCG  
 CGTCGGGCCACCACTGCTCATCGAGCATCTGGGGCTGGCGGGGGCGGCGCAGCAAGACTGCGCCTG  
 TCGTGGGCTGCGGCTGGGTGCGCGGCCGCTGCGGGCTCCCGCCGGGGTCCCACCGCGCTGCCTGCCT  
 ACCCGGCGGCCGAGCCGGGGCCGCTGTGGCTGCAGGGCGAGCCTCGTCACTTCTGCTGCCTGGACTTCAG  
 CCTGGAGGAGCTGCAGGGCGAGCCGGGCTGGCGGCTGAACCGCAAGCCATCGAGTCTACGCTGGTGGCC  
 TGCTTCATGACCTGGTCATCGTGGTGTGGAGCGTGGCCGCCCTCATCTGGCCGGTGGCCATCATTGCCG  
 GCTTCTGCCAACGGCATGGAGCAGCGCCGACCACCGCCGGCGCCCCGCGGGCGCTCCTGCAGCCGT  
 ACCCGCGGGGACTACCGCCGCTGCCCGCCGCTGCTGCAGCTGCTGCTGCTGCCCGGCTGTCAAC  
 TCGGGGTGGCTCCTAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG216124 representing NM\_001002267  
 Red=Cloning site Green=Tags(s)

MLPLLAALLAAACQLPPAHGGATDAPGLAGTPPNASANASFTNEHSTPRLLASAASAPPERAGPEEAPAA  
 PCNISVQRQMLSSLLVRWGRPRGLQCDLLLFTNAHGRAFFAAAFHRVGPPLLIEHLGLAAGGAQQDLRL  
 CVGCGWVRGRLRAPAGAPTALPAYPAEPLWLGEPHFCCCLDFSLLEELQGEPLNRKPIESTLVA  
 CFMTLVIVVWSVAALIWPVPIIAGFLPNGMEQRRRTTAGAPAAAPAAVPAGTTAAAAAAAAAAAAAAAAAVT  
 SGVAPK

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001002267

**ORF Size:** 858 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001002267.2](#), [NP\\_001002267.2](#)


**RefSeq Size:** 1744 bp

**RefSeq ORF:** 861 bp

**Locus ID:** 72309

**Cytogenetics:** 9 F4

### Product images:



Circular map for MG216124