

## Product datasheet for **MG216178**

### Ssr4 (NM\_009279) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ssr4 (NM\_009279) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ssr4  
**Synonyms:** AW743063; C78062; Trap  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG216178 representing NM\_009279  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCGGCGATGGCATCTTTGGCGCCCTGGCACTACTCCTGCTGTCTGGCCTATCTTGCTGCTCAGCAG  
 AGGCCTGCCTGGAACACAGATCACCCCTTCTTACTATAACAATTCAGATGCCGTCATTTCTACAGAGAC  
 TGTATTCATCGTGGAGATTTCACTGACCTGCAAGAACAGGGTCCAGAACATGGCTCTTTATGCCGACGTT  
 AGTGGAAAACAATTTCTGTAACCCGCGCCAGGATGTGGGTCGATATCAGGTGCTCTGGAGCCTGGAGC  
 ACAAGAGCGCTCATGCGGGCACCTATGAGGTCAGATTCTTTGATGAAGAGTCCTATAGCCTCCTAAGGAA  
 GGCTCAGAGAAAATATGAGGACGTTTCTATCATTCCACCCTGTTACAGTCAGTGTGGACCATCGGGGT  
 ACCTGGAATGGCCTTGGGCTCTCCACAGAGGTGCTGGCTGCAGTAATTGGCATAGTGATCTACTACCTAG  
 CTTTCAGTGCAAAGAGCCACATCCAGGCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG216178 representing NM\_009279  
 Red=Cloning site Green=Tags(s)

MAAMASFGALALLLSGLSCCSAEACLEPQITPSYYTTSDAVISTETVFIVEISLTCKNRVQNMALYADV  
 SGKQFPVTRGQDVGRYQVSWLSLEHKS AHAGTYEVRFFDEESYLLRKAQRNEDVSIIPPLFTVSVDRHG  
 TWNGPWVSTEVLA AVIGIVIIYLA FSAKSHIQA

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_009279

**ORF Size:** 516 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_009279.3](#), [NP\\_033305.1](#)

**RefSeq Size:** 703 bp

**RefSeq ORF:** 519 bp


**Locus ID:** 20832

**UniProt ID:** [Q62186](#)

Cytogenetics: X 37.41 cM

Gene Summary: TRAP proteins are part of a complex whose function is to bind calcium to the ER membrane and thereby regulate the retention of ER resident proteins.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MG216178