

Product datasheet for **MG216218**

Samd3 (NM_001115154) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Samd3
Synonyms:	Gm623
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >MG216218 representing NM_001115154
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAACGTGGTCAGTTGATCAGGTCTGCAAATGGTTGGTGGAAAAAACTTAGGTGAGCTAGTCCCA
 GGTTTCAAGAGGAAGAAGTCAGTGGGGCCACACTCCTGGCCCTTAATGATCGGATGGTACAGCAACTGGT
 AAAGAAAAATGGACACCAGGCTGTTCTGATGGACTTCATTAAGAAATACAAGCAGGGCAATCAAGAGCTG
 AAGCCTACGGGAGGCCCGCAGATACATCCACGCTGACCCCGGCACAAGCGGCCCTGAGCACGAACAAA
 ACCCTAGTCCCACCAGCCACGGGGACCAGACATCTTTGTATCCAGCCGTCCTTGACAATAGACTCATTGA
 TCAAAGAGTATTGAAGCAGAGAAGAAATGTAACAATGTTCTAGCAAGGCATAAAGCATTACAGTGGACT
 AAGTCCTACATCTTACCGAGTTCCCTATGATGTCAAGTGCATGTTGGTGGAGCAGAAGCGCCCGACC
 ACAGCATGCGGATAAAGATCATTGAGTTTCTCCAGGCAGATAGACTAAGTACCTGGAAGGCTCGCTGTA
 CCCCACCACCAGCAGTACAACGACGTGGTCAATGCTCTGCTGCAGGCCACCCTTTCTGGATGAGGAC
 GGCTGTGGCTTTTTTTATGGAACAGCCCTCAAAGATCGCTTTAAATATATTGAAAGACCATAGAAG
 ATGATGAACAAGTAATGAGGAATAAGTGTAAATTTGGACACCGAAGAGGCCAGACAAGGAAGTCTCTTG
 TGACATACAATCGAATGAAATTAATAATTTGTGCAGATCAAAGAAGAATCTGCTCATTTAGACTCGGAGTG
 GATGAACACATCAGCTGGTCCAGCAAGAGTATATGAAAACAGAAAGGGACTGGAGAGAAGTTGACAAGA
 GGATGAGCCAGACTTTGGAAATAAGAAGAAGATGATTGGCGGCCAAACACCTCTCAAGGACATCTTAA
 GATGTTTCTTCTGAAAGTGCCTTACCAGATGTTCCAGAGAAGTTCAAATTTTACAAAAACAGATATT
 TATAAGAAAACAAGGCATATCCTGGAATCTACTCGGAAAAATACTGACTGCCTTTTCAGTGTGGACA
 ATCCAATCAATACTGCCCTGCAAGAGAAGATGAAACACTACACTGACGAAGGGTGTGAAGTACATGAA
 GATGACTGCCACGTGTTTACTCCTCCACATGTTTTGGGGATGAGCCAGCCTTTTTGTCGTTGTGAAT
 GGAAGGTGCACGTGCCACGCTGTGTTGGAAGTAAAAACCTTTCCACATCAATGGCTGTGAGTTTT
 CTCTGTATTTGAACAAAGAGAAGCTCACAAGGTGGATGACTGTGTACAGCTCGGCTGCTTGGTGTG
 TGCTTTCCGTGTGTTGGGATCGAGTGTCCAAGAAGACTGTCCAAACCTTCACTTCTAGAGACATTG
 ATCTTTGACATGCAGAGTCCCAAGTTCCCTCTTTGAAAGAAAAGGAAATCAGATCTCAACCTCCAATCA
 CT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >MG216218 representing NM_001115154
 Red=Cloning site Green=Tags(s)

METWSVDQVCKWLVEKNL GELVPRFQEEVSGATLLALNDRMVQQLVKKIGHQAVLMDFIKYYKQGNQEL
 KPTGGPADTSTLTPAQAAPEHEQNPSPTSHGDQTSLYPAVLNRLIDQRLVKQRRNVKHLARHKALQWT
 KSYILPEFPYDVKMLVEQKRPDHSRIRIIEFLQADMTKYLEGSLYPTTQYNDVNVNALLQHPFLDED
 GCGFFLWKRALKDRFKYIRRP IEDDEQVMRNKCKFGHRRGQTRKSLADISNEIKIVQIKEESAHL DSEV
 DEHISWFQEQEYMKTERDWREVDKRMSTLEIRRKMI GGQTP LKDI LKMF PFLKCPYQMFREVQILTKTDI
 YKTRHILESYSENILTAFSVLDNPINTALQEKMKHYTDEGLKYMKMTATCLLLPHVFGDEPSLFV VVN
 GKVVHSTPVLVKNPFHINGCEFSLYLNKEKLT KVDDCVTALAALVSAFRVFGIECPRRLSQTFNFLETL
 IFDMQSPQFP SLKEKEIRSQPPIT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001115154

ORF Size: 1542 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001115154.1](#), [NP_001108626.1](#)

RefSeq Size: 1892 bp

RefSeq ORF: 1545 bp
 Locus ID: 268288
 UniProt ID: [Q8C4H2](#)
 Cytogenetics: 10 A4

Product images:



Circular map for MG216218