

Product datasheet for **MG216225**

Ccdc50 (NM_026202) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc50 (NM_026202) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ccdc50
Synonyms:	C3orf6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG216225 representing NM_026202 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCGACGTGAGTGTAGATCAGTCCAAGTTGCCGGCGTGAAGGAAGTGTGTAGAGATTTTGCCGTCC
TGGAGGACCACACCCTGGCCATAGCCTGCAGGAACAAGAGATTGAGCATCATTTGGCATCCAACATTCA
GCGGAACCGTCTGGTACAACATGATCTGCAGGTTGCTAAGCAGCTCCAAGAGGAAGACCTCAAAGCCCAA
GCTCAGCTCCAGAAGCGCTACAAAGCCCTGAACAACATGATTGTGAAATTGCTCAGGAAATCCAGGAGA
AGCTAACCATTTGAGGCTGAGAGACGACGCATTCAGGAAAAGAAGGATGAGGACATAGCACGTCTTTTGCA
AGAGAAGGAGCTACAGGAGGAGAAAAGGAGGAAGAAACACTCCAGAGTTTTCTGGGGCAGTGTTTT
GGAGATAACTACTATCATGAGGATGGAGGAATGAAACCAAGAGGAATAAAGAAGCTGTCTACTCCAG
CACGAGCAAGCCACAGGGACCAGGAGTGGTATGATGCTGAAATTGCCAGGAAATTGCAAGAAGAAGAACT
TTTGGCTACTCATGTGGACATGAGAGCAGCTCAGGTTGCCAGGATGAGGAAATTGCTCGACTTCTAATG
GCTGAAGAAAAAAGCTTACAAGAAAGCCAAAGAGCGAGAAAAGTCATCTTTGGACAAAAGGAAACATG
ACCCTGAATGCAAGTTAAAAGCAAAGTCAGCCCACTCAAAGTCAAAGAGGGTGATGAAGCACACCCGCTC
CAAGATTGACAGGCCATCAAGACCACCACCTACCATGATGGGCTTGAGGACACAGATCCCACCCAT
TTTACAAACCAGCACAGTACAACATGGCATCTTCAAAGTCAGAGTCCTCACAGAAAGGCTTCCATAACA
AGCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG216225 representing NM_026202
 Red=Cloning site Green=Tags(s)

MADVSDVQSKLPGVKEVCRDFAVLEDHTLAHSLQEQEIEHHLASNIQRNRLVQHDLQVAKQLQEEDLKAQ
 AQLQKRYKALEQHDCEIAQEIQEKL TIEAERRRIQEKDEDIARLLQEKLQEEKRRKKHTPEFSGGSVF
 GDNYYHEDGGMKPRGIKEAVSTPARASHRDQEWYDAE IARKLQEEELLATHVDMRAAQVAQDEE IARLLM
 AEEKKAYKKAKEREKSSLDKRKHDPCKLKAKSAHKSKEGDEAHRSKIDRPSRPPPTMMGLETDPTH
 FTNQHSTTWHLPKSESSQKGFHNKQ

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_026202

ORF Size: 915 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_026202.4 , NP_080478.2
RefSeq Size:	6981 bp
RefSeq ORF:	918 bp
Locus ID:	67501
UniProt ID:	Q810U5
Cytogenetics:	16 18.98 cM
Gene Summary:	Involved in EGFR signaling.[UniProtKB/Swiss-Prot Function]