

Product datasheet for **MG216237**

Dph3 (NM_172254) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dph3 (NM_172254) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Dph3
Synonyms: DELGIP1; DelgipP1; Desr1; KTI11; Zcsl2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG216237 representing NM_172254
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGTGTTTCACGACGAGGTGGAGATCGAGGACTTTCAATATGACGAGGACTCGGAGACATATTTCT
ACCCTTGCCCTGTGGGATAACTTTGCCATACCAAGGAAGATTTGGAAAATGGAGAAGATGTGCCAC
GTGTCCTAGCTGCTCACTCATTATAAAGTGATTTATGACAAAGATCAGTTCATGTGTGGAGAAACAGTC
CCAGCACCTTCAACCAACAAGGAGTTAGTAAATGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG216237 representing NM_172254
Red=Cloning site **Green**=Tags(s)
MAVFHDEVEIEDFYDEDESETYFPCPCGDNFAITKEDLENGEDVATCPCSLIIKVIYDKDQFMCGETV
PAPSTNKELVKC

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_172254

ORF Size: 246 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172254.4
RefSeq Size:	3016 bp
RefSeq ORF:	249 bp
Locus ID:	105638
UniProt ID:	Q8K0W9
Cytogenetics:	14 B
Gene Summary:	Essential for the first step in the synthesis of diphthamide, a post-translational modification of histidine which occurs in elongation factor 2.[UniProtKB/Swiss-Prot Function]