

## Product datasheet for **MG216277**

### Ifne (NM\_177348) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Tag:** TurboGFP  
**Symbol:** Ifne  
**Synonyms:** Ifn-tau-1; Ifne1; Ifnt1; Infel  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG216277 representing NM\_177348  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGTTCACAGACAGCTCCCTGAAACGGTGTGCTGCTCTTGGTTTCTTCCACTATCTTCTCCCTAGAAC  
CGAAACGGATTCCCTTCCAATTGTGGATGAACAGAGAAAGCCTACAACACTGAAACCTTTGCCAAGCTC  
GTCAGTCCAGCAGTGTCTAGCACACAGGAAGAATTTCTGCTTCTCAGCAGCCTGTGAGTCTCACCAG  
TACCAAGAGGGACAGGTGCTGGCTGTTGTGCACGAGATCCTTCAGCAGATCTTCACGCTCCTCCAGACAC  
ATGGGACTATGGGCATTTGGGAGGAAAACCATATAGAAAAAGTCTTAGCTGCGCTTACCCGGCAGCTGGA  
ATACGTGGAGTCACTGGGTGGACTGAACGACGCGCAGAAGAGTGGGGGCTCGAGTGCCGAGAACCTTAGG  
TTACAGATTAAGCATACTTCAGGAGGATCCACGATTACTTGAAAACCAAGGTACAGCAGCTGTGCCT  
GGATCATTGTCCAGACAGAAATCCACCGCTGTATGTTCTTTGTGTTTCAGTTTCAACATGGCTGAGCAG  
ACAAGACCCAGACCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG216277 representing NM\_177348  
Red=Cloning site Green=Tags(s)

MVHRQLPETVLLLLVSSTIFSLEPKRIPFQLWMNRESLQLLKLPLSSSVQQCLAHRKNFLLPQQPVSPHQ  
YQEGQVLAVVHEILQQIFLLQTHGTMGIWEENHIEKVLAAALHRQLEYVESLGLNAAQKSGSSAQNLRL  
LQIKAYFRRIHIDYLENQRYSSCAWIIVQTEIHRMFFVFRFTTWLSRQDPDP

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NM\_177348

**ORF Size:** 576 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_177348.2](#), [NP\\_796322.1](#)

**RefSeq Size:** 664 bp

RefSeq ORF: 579 bp


Locus ID: 230405

UniProt ID: [Q80ZF2](#)

Cytogenetics: 4 42.04 cM

**Gene Summary:** Type I interferon required for maintaining basal levels of IFN-regulated genes, including 2'-5'-oligoadenylate synthetase, IRF7 and ISG15, in the female reproductive tract. Directly mediates protection against viral, including HSV-2, and bacterial, including Chlamydia muridarum, genital infections.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MG216277