

Product datasheet for **MG216332**

Nme5 (NM_080637) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nme5 (NM_080637) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Nme5
Synonyms: 1700019D05Rik; Nm23-M5
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG216332 representing NM_080637
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGTATCAATGCCCTGCCTCAGATATATGTAGAAAAACCCCTAGCCCTTATCAAGCCAGATGTTG
TTGACAAAGAAGAGGAGATACAGGATATTATTCTGGGATCTGGATTCACCATTATTCAGAGACGGAACT
ACACCTGAGCCCCGAGCACTGTAGTAACTTTTATGTGGAGCAGTATGGGAAAATGTTTTCCAACTTA
ACAGCTTATATGAGCTCTGGACCTCTTGTGCTATGATATTAGCTAGACATAAGGCCATCTCCTACTGGA
AAGAACTTATGGACCAAGTAACAGTTTAGTAGCCAAGGAGACACACCCGGACAGCTTAAGGGCGATATA
CGGTACAGACGAGCTCAGGAACGCACTTCATGGGAGCAATGACTTCGCTGCCTCAGAGCGAGAGATCAGG
TTCATGTTTCCAGCCGTGATTATTGAACCCATTCCAATTGGACAAGCCGCTAAGGACTACATAAACTTGT
ACGTAGCGCCAACCCTACTTCAAGGACTCACAGAGCTTTGTAAGGAAAAGCCACCAGACCCTTACCTCTG
GCTCGCTGATTGGCTGATGAAAAATAATCCAACAACCTAAACTTTGTCATTTCCAGTGACAGAAGAG
CCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG216332 representing NM_080637
 Red=Cloning site Green=Tags(s)

MEVSMPLPQIYVEKTLALIKPDVVDKEEEIQDIIILGSGFTIIQRRKLHLSPEHCSNFYVEQYGMFFPNL
 TAYMSSGPLVAMILARHKAIISYWKELMGPSNSLVAKETHPDSLRAIYGDEL RNALHGSNDF AASEREIR
 FMFPAVIIIEPIPIGQAAKDYINLYVAPTL LQGL TELCKEKPPDPYLWLADWLMKNNPNKPKLCHFPVTEE
 P


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_080637

ORF Size: 633 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_080637.4
RefSeq Size:	829 bp
RefSeq ORF:	636 bp
Locus ID:	75533
UniProt ID:	Q99MH5
Cytogenetics:	18 B1
Gene Summary:	Does not seem to have NDK kinase activity. Confers protection from cell death by Bax and alters the cellular levels of several antioxidant enzymes including Gpx5. May play a role in spermiogenesis by increasing the ability of late-stage spermatids to eliminate reactive oxygen species.[UniProtKB/Swiss-Prot Function]