

## Product datasheet for **MG216387**

### **Mcemp1 (NM\_026985) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Tag:** TurboGFP  
**Symbol:** Mcemp1  
**Synonyms:** 1810033B17Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG216387 representing NM\_026985  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGCATGCATCAGCCTCCAGGATAAGAACCGGAGGAAGCCAGGTCATGATGAAGGTGCTACAATCCTG  
ACTACGAGAATATAACCTTGGCCTTCAGAAACAAGGACCAACTCAAACCTCAGCCAATCAACACCCACAAA  
ACAAGCCAAGTTCAAGACATCCCTGGACCCAGCTGAGTCCCGCCTTGGTTGTACAGAACCATTATGATG  
TTGTATGTTCTCCTTGCCTCGCTTTTTATCCTGCATCGTCCTCTGCTTTGGTCTTGGTAAAAAATT  
CTGAGATGTCCAAGGAGCTGTGGACCTTGAAGCAGAGCTTTCGAATGTTTCAGACACGGTGTGGAATAT  
CCGGGAGCTCCAGAATCAGCAAACGAGGATTTGGGAAGCTGCCAGGGGACATCAAGGAGGTCAAGAAG  
ACCCTTGGCACAGTCATGAGTAGCATCCAGACTGAAAACGACCGGCTGAAGACTGTGCCGGCAGATATAA  
CCCAATCAAGAAAACCTTGAGGCGCTAGAAAAGAAGGCACAGCCTCAGCCAGTACA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG216387 representing NM\_026985  
Red=Cloning site Green=Tags(s)

MHASASQDKNRRKPGHDEGAHNPDIENITLAFRNKDQLKLSQSTPTKQAKFKTSLDPAESPWLVRTIMM  
LYVLLALVFLSCIIVLSALVLVKNSEMSKELWTLKAE LSNVSDTVWNIRELQNNQTRIWEAAQGDIKEVKK  
TLGTVMSSIQTGNDR LKTVPADITQIKKTLEALEKKAQPQPST

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NM\_026985

**ORF Size:** 549 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_026985.1](#), [NP\\_081261.1](#)

RefSeq Size: 1132 bp  
 RefSeq ORF: 552 bp  
 Locus ID: 69189  
 UniProt ID: [Q9D8U6](#)  
 Cytogenetics: 8 1.92 cM

**Product images:**



Circular map for MG216387