

Product datasheet for **MG216469**

Tifab (NM_001168615) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tifab (NM_001168615) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Tifab
Synonyms: MGC32345; MGC37193
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG216469 representing NM_001168615
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGAGGCCCTCACAGTCTGCAGGTGAGCCTGTACCACCCTACACAGGGCCCGGTGCGCTTTGCC
 ACGTCCCACAGCAGCTGCAACATGACGCCAGCCGGCTGTTGGTTGGGCGAGGACAGAACACCCATCTCCA
 GCTGCAGTGCCACAGTTGTCTCGCTACCATTTGTCCTGGAACCCTACCTGGAGAAAGGCAGCAGTCTG
 CTGGCCTTCTGCCTCAAGGTGCTGACCCGCAAGAGCTGTGTGGGTCAACGGGCTGCCACTGAGGTACC
 TGGAGCAGTTCCCTTGGGCACCATCAACAGAATCTCCTTCTCTGGCATCCAGATGCTAGTCCGAAAGA
 GGGGGGAGCTCCCTTGAGACCTTTGTCTGCTACTTCCATCTCAGCCCTCACCCCTGATTACAGACCC
 AAGGCTCAGGAGACAGATGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG216469 representing NM_001168615
 Red=Cloning site Green=Tags(s)

MERPLTVLQVSLYHPTQGPVAFVAFVHPVQQLQHDASRLLVGRGQNTLQLQLPQLSRYHLSLEPYLEKGSLL
 LAFCLKVLTRKSCVWVNGPLRYLEQVPLGTINRISFSGIQMLVRKEGGASLETFCYFHLSPSPLIYRP
 KAQETDE

TRTRPLE - GFP Tag - V


Restriction Sites: Sgfl-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_001168615

ORF Size: 441 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001168615.1, NP_001162086.1</u>
RefSeq Size:	3498 bp
RefSeq ORF:	444 bp
Locus ID:	212937
UniProt ID:	<u>Q8JZM6</u>
Cytogenetics:	13 B1
Gene Summary:	Inhibits TIFA-mediated TRAF6 activation possibly by inducing a conformational change in TIFA. [UniProtKB/Swiss-Prot Function]