

## Product datasheet for **MG216476**

### **Olf1102 (NM\_207154) Mouse Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Olf1102
Synonyms:	MOR179-4
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG216476 representing NM\_207154, **codon optimized**.  
Due to the complexity of NM\_207154, the ORF clone is codon optimized for mammalian Expression.  
The nucleotide sequence differs from the reference sequence, yet the amino acid sequence remains identical.

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCCAAGAACTCCTTCTTACACTAACACAAAAACAACACTCAGGTGAATAACGTGACAGAGATCACAGTGT  
TCATCTGCTCGGATTTACAGATGACGTTGACATGAATATCTTTCTGTTTCATCCTTTTTTTGGCAATCTA  
CGTGGTGACTTTGATCGGGAATTTGGGCCCTCGTCGTGCTGGTGATAGAAGACTCTCGATTGCATAACCCC  
ATGTAATACTTTCTCACAGTCTGTCTCCCTCGACGCATGTTTCTCATCTGTGTTGACACCAAAAATGT  
TGGCTAATTTTCTGAGCAAGAACAAGATATCTCCTTTGCCGGCTGCGCCACCCAGATGTTGTTGTTTGT  
CACATTCGGCACCACTGAATGTTTCTGCTCGCGGCAATGGCATAACGACAGATACCTCGCTATCTACAGC  
CCACTTCTGTATGCCGTCGGATGTCACCCAGAGTGTATGTTCCGTTGATCATCGCCAGCTACACCGGTG  
GAATCCTTCACGCCACTATCCACACTGTGGCCACCTTTCCCTGAGTTTCTGCGGTAGCAACGAGATAAG  
GCACGTTTTCTGCGATATCCCTCCCTGTTGGCGCTGTCTTGAGCGATACCCATCTGAACCAAGCTGCTC  
CTGTTCTACTGCGCCGGTTCCATTAAGCTGATTACCATACTGATAGTCTGTTAGCTACGGGTTTGTCC  
TTCTCGCAATCCTCAAATAAATAGCGCTGAAGGCAGAAGAAAGATATTCTCCACTTGGGGGCCACCT  
GACTGGTGTGCAATCTTTCACGGCACCATTCTGTTTCATGTACGTGCGGCCCTTCAAACCTACACCTG  
GAGCAGGACATGGTGGTTTCTACATTCTACACCATCGTGATTCCAATGCTTAATCCTATCATATATAGCC  
TGAGAAAACAAAGACGTGAAGGAGGCGATGAGAAAATTGCTCAAGAGAAAACCTTGCCACGAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



**Protein Sequence:** >MG216476 representing NM\_207154  
Red=Cloning site Green=Tags(s)

MPRTPSYTNKTTQVNNVTEITVIFILLGFTDDVDMNIFLFLFLAIYVVTLIGNLGLVVLVIEDSRLHNP  
 MYYFLTVLSSLDACFSSVLTPKMLANFLSKNKSISFAGCATQMLLFVFTGTTECFLLAAMAYDRYLAIYS  
 PLLYAVRMSPRVYVPLIIASYTGGILHATIHVATFSLSFCGSNEIRHVFCDIPLLALSCSDTHLNQLL  
 LFYCAGSIKLITILIVLSYGFVLLAILKINSAEGRRKIFSTCGAHLTGVSIFHGTLIFMYVRPSSNYTL  
 EQDMVVSTFYTIIVPMLNPIIYSLRNKDVKEAMRKLKRKLVHE

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**ACCN:** NM\_207154

**ORF Size:** 972 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<u>NM_207154.2, NP_997037.2</u>
<b>RefSeq Size:</b>	1063 bp
<b>RefSeq ORF:</b>	975 bp
<b>Locus ID:</b>	228228
<b>UniProt ID:</b>	<u>Q8VES2</u>
<b>Cytogenetics:</b>	2 D
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG216476