

Product datasheet for **MG216560**

Puf60 (NM_001164600) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Puf60 (NM_001164600) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Puf60
Synonyms:	2410104I19Rik; 2810454F19Rik; SIAHBP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG216560 representing NM_001164600
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAAAATGGACAAAGCACAGGCACCAAGCTGGGGCTGCCTCCCCTGACGCCCGAGCAGCAGGAGGCC
 TCCAGAAGGCCAAGAAATATGCAATGGAGCAGAGCATCAAGAGTGTGCTGGTGAACAGACCATCGCGCA
 CCAGCAGCAGCAGCTCACCAACCTGCAGATGGCAGCTCAACGGCAGCGGGCACTGGCTATCATGTCCGG
 GTCTATGTGGGCTCCATCTACTATGAGCTGGGAGAAGACACTATTCGCCAGGCCTTTGCTCCCTTTGGCC
 CCATCAAGAGCATTGATATGTCTGGGACTCCGTTACCATGAAGCATAAGGGCTTTGCCTTCGTGGAGTA
 CGAGGTCCAGAAGCTGCGCAGCTCGCTTTGGAGCAGATGAACCTGTGTATGCTTGGGGCAGGAACATC
 AAGGTGGGAAGACCCAGCAACATTGGACAAGCCCAACCCATCATAGATCAGCTGGCTGAGGAAGCTAGGG
 CTTTCAACCGAATATATGTGCCCTCTGTACATCAAGACCTGTCTGATGATGACATCAAGAGTGTATTTGA
 AGCCTTTGGCAAGATCAAGTCTTGTACGCTGGCCCGGGACCCACAACCTGGCAAGCACAAGGGCTATGGT
 TTTATCGAATATGAGAAGGCCAGTCGTCCCAGGATGCTGTGTCTCCATGAACCTCTTTGATCTGGGTG
 GCCAGTACTTGAGGGTGGCAAGGCCGTACACCCCCATGCCCTGCTAACACCTGCCACGCCTGGAGG
 TCTCCCGCTGCTGCTGTGTGGCCGAGCTGCAGCCACGGCCAAGATTACAGCTCAGGAAGCAGTGGCT
 GGAGCAGCAGTGTGGTACTTTAGCCACACCAGGACTAGTGTCCCAGCACTGACACTGGCCAGCCCT
 TAGGGGCCCTACCCAGGCTGTATGGCTGCCAGGCCCTGGAGTCATCACAGGTGTGACACCAGCCCG
 CCCTCCTATTCCGGTACCATCCCCTCTGTGGGAGTAGTGAACCCATTCTCGCCAGCCCGCCAACGTTG
 GGTCTGTTGGAGCCCAAGAAGGAGAAGGAAGAGGAGGAGCTGTTCTCTGAGTCAGAACGGCCAGAGATGT
 TGAGTGAGCAGGAGCACATGAGCATCTCGGGCAGCAGTGTCTGCCACATGGTCATGCAAAAAGCTACTCAG
 AAAACAGGAGTCTACAGTGTGTTCTCCGAAACATGGTGGACCCCAAGGACATCGATGACGACCTGGAG
 GGAGAGGTGACAGAGGAGTGTGGCAAAATTTGGTCTGTGAACCGGTATTATCTATCAAGAAAAACAGG
 GCGAAGAGGAAGATGCAGAAATCATCGTCAAGATTTTTGTGGAATTCTCCATGGCCTCAGAGACTCACAA
 GGCCATCCAGGCCCTCAATGGGCGCTGTTTGGTGGTCCGAAGTGGTGGCTGAAGTGTATGACCAGGAG
 CGTTTTGATAACAGCGACCTCTCCGCG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG216560 representing NM_001164600
 Red=Cloning site Green=Tags(s)

MENGQSTGTKLGLPPLTPEQQEALQKAKKYAMEQSIKSVLVKQTIHQQQQLTNLQMAAQRQALAIMCR
 VYVGSIIYELGEDTIRQAFAPFGPIKSIDMSWDSVTMKHKGFVVEYEVPEAAQLALEQMNSVMLGGRNI
 KVGSRPSNIGQAQPIIDQLAEERAFNRIVVASVHQDLSDDDIKSVFEAFGKIKSCTLRDPTTGKHKGYG
 FIEYEKAQSSQDAVSSMNLFDLGGQYLRVGKAVTPPMPLLPATPGGLPPAAAVAAAAATAKITAQEAVA
 GAAVLGTLATPGLVSPALTLAQPLGALPQAVMAAQAPGVIITGVTPARPPIPVTIPSVGVVNPILASPPTL
 GLLEPKKEKEEEEELFPESERPEMLSEQEHMSISGSSARHMVMQKLLRKQESTVMVLRNMVDPKDIDDDLE
 GEVTEECGKFGAVNRVIIYQEKQGEEDAEIIVKIFVEFSMASETHKAIQALNGRWFGRKVVVAEYVDQE
 RFDNSDLA

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001164600

ORF Size: 1497 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164600.1](#), [NP_001158072.1](#)

RefSeq Size: 1848 bp

RefSeq ORF: 1500 bp

Locus ID: 67959


UniProt ID: [Q3UEB3](#)

Cytogenetics: 15 D3

Gene Summary:

DNA- and RNA-binding protein, involved in several nuclear processes such as pre-mRNA splicing, apoptosis and transcription regulation. In association with FUBP1 regulates MYC transcription at the P2 promoter through the core-TFIID basal transcription factor. Acts as a transcriptional repressor through the core-TFIID basal transcription factor. Represses FUBP1-induced transcriptional activation but not basal transcription. Decreases ERCC3 helicase activity. Is also involved in pre-mRNA splicing. Promotes splicing of an intron with weak 3'-splice site and pyrimidine tract in a cooperative manner with U2AF2. Involved in apoptosis induction when overexpressed in HeLa cells. Modulates alternative splicing of several mRNAs. Binds to relaxed DNA of active promoter regions. Binds to the pyrimidine tract and 3'-splice site regions of pre-mRNA; binding is enhanced in presence of U2AF2. Binds to Y5 RNA in association with TROVE2. Binds to poly(U) RNA (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG216560