

Product datasheet for **MG216632**

Xlr3b (NM_001081643) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: Xlr3b
Synonyms: EG574437; Xlr3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG216632 representing NM_001081643
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCAAGCAGAAAAAGGAAGGCCACTGACACTGCTGGCAGGCACTCAAGGATGGATCCAAATCTCTCCT
CAGATGACAGCCAGAACCAGGTGCAGTCGAGCAGCTAACAGAGAAGTCTTGATGCTGGTAGGGAGGA
CATTATTTCTCAGGAACAGAGAGGCCAACAGGCCAGGAAGGAAAAACAGGATTTGGTCCAAGAATTCGAA
GAGCCACCTAACAAAGTTCTTCAGGAAAAACAGAGAAAGATTCTCAAGGATCATGACCTCATCTTTTCAGTG
CCATGGAGGTCAAATTAAGGATGTTCTGAAAACCCACTGTGAAGAGAGAGGCAGAAACTTTGTCAGGACTA
TTCTCTTCAGTTTACAAATTTGAATAGGAAGTTAACTTCGGATGCATACAACTCAAGAAACATGCCGAA
ACACTCTAATATGTTTATGGAGCAACAGAAGTTTATTCATGAAAGTCTCACTCTCAGAAGAACAGAA
TGAAGGAATTTAAATCACTGTGTGAAAAACTTGGAGAACTGGAGGTACTGAGGGATTCTCGGGGAAA
TTCCATCGCTGAAGAGCTGAGACGTCTCATAGCCACCTTGGAAATCAAATCTGATGCTGCATAACCAG
CAAAAGACTGCTGCTCCTCCACAGTCTCTCTGGATGTGTTATTCTCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >MG216632 representing NM_001081643
Red=Cloning site Green=Tags(s)

MSSRKRKATDTAGRHSRMDPNLSSDDSQNPGAVAAANREVL DAGREDI I SSGTERQQARKEKQDLVQEF E
PPNKVLQENRERFSRIMTSSF SAMEVKIKDVLKTHCEERQKLCQDYSLQFTNLNRKLTSDAYLKKHAE
T LSNMFMEQQKFIHESL TLQKNRMKEFKSLCEKYLEKLEVL RDSRGN SIAEELRRLIATLEIKLLMLHNQ
QKTAAPPQSLLDVLF S

TRTRPLE - GFP Tag - V



RefSeq Size:	1636 bp
RefSeq ORF:	681 bp
Locus ID:	574437
UniProt ID:	Q6P205
Cytogenetics:	X A7.3

Product images:

Circular map for MG216632