

Product datasheet for **MG216646**

Prss33 (NM_001081399) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prss33 (NM_001081399) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Prss33
Synonyms: Eos; mT6; tryptase-6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG216646 representing NM_001081399
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGGGGTGCTCCACCTCCAGATCCTGCTCCTGTTGGTGTTAGGAACCAGGATGCAGGAGTGTGCAG
 CTTGTGGGCAGCCCCGCATGTCCAGCCGGATCGTTGGGGCCGGGATGCTCAGGATGGAGAGTGGCCGTG
 GCAAACAAGCATTAGCACCCTGGAGCCCATGTGTGCGGAGGGTCACTCATTGCACCCAGTGGGTGTTG
 ACGGCAGGTCACTGTTTCCCTAGACGGGTATGGCCGTGAGTACAGTGTGCTCCTGGGAGCACTGAGTC
 TGGACGTGAGTCACTCCACGAGCTCTGGTTCCTGTGCTGCGGGTGTGCTGCCTCCTGACTACTCTGA
 GGACGAGGCCCGTGGTGACCTGGCACTGCTGCAGCTGCGTCACCCAGTGTCCCTGAGCACCCGCATCCAG
 CCTGTCTGCCTGCCTGCACCAGGTCCCATCCACCACCAGGGTCTCCATGCTGGGTCACTGGCTGGGGTA
 GCCTCAGCCCAGGAGTACCACTCCAAAGGGGAGACCCTTGAAGGAGTGAAGGTGCCGCTGCTGGACTC
 CCGAGCCTGTGACCGCCTTACCACGTGGGTGCCAACGTGCCAGGGTGAACGCATAGTGTACCAGGG
 AACTTGTGTGCTGGTACCGCAGAGGCCACAAGGATGCCTGCCAGGGTACTCTGGAGGACCCTGACCT
 GTATGGAGTCTGGCCACTGGTCTGGTGGTGTGGTGAAGTGGGCAAGGGCTGTGCCCTGCCAACCC
 TCCAGGTGTCTATACCAATGTGGCAAGTATAGTCCCTGGATCCAGGCTCGCCTGAGCCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG216646 representing NM_001081399
 Red=Cloning site Green=Tags(s)

MRGASHLQILLLLVLGTRMQECAACGQPRMSSRIVGGRAQDGEWPWQTSIQHRGAHVCGGSLIAPQWVL
 TAGHCFPRRVWPSEYSVLLGALSLDVRSSHELLVPVLRVLLPPDYSEDEARGDLALLQLRHPVSLSTRIQ
 PVCLPAPGSHPPPSPCWVTGWGSLSPGVPLPKGRPLQGVRVPLLDLRACDRLYHVGANVPQGERIVLPG
 NLCAGYRRGHKDACQGDSSGGLTCMESGHVVLGVVSWGKGCALPNRPGVYTNVAKYSPWIQARLSL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001081399

ORF Size: 831 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001081399.2 , NP_001074868.2
RefSeq Size:	1384 bp
RefSeq ORF:	834 bp
Locus ID:	353130
UniProt ID:	Q80WM7
Cytogenetics:	17 12.08 cM
Gene Summary:	Serine protease that has amidolytic activity, cleaving its substrates before Arg residues. [UniProtKB/Swiss-Prot Function]