

Product datasheet for **MG216709**

Steap2 (NM_028734) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Steap2 (NM_028734) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Steap2
Synonyms:	4921538B17Rik; AI930049; AW045895; IPCA-1; IPCA1; PCANAP1; STAMP1; STMP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG216709 representing NM_028734
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAATCCATCTCTATGATGGGAAGCCCTAAGAGCCTGGAGACGTTTTTGCCTAATGGCATAAATGGTA
 TCAAAGACGCAAGGCAGGTCACCGTGGGGTGATAGGAAGTGGGATTTTGCCAAGTCTCTGACCATTCCG
 GCTTATTAGATGTGGCTACCACGTGGTCATAGGAAGCAGAAATCCCAAGTTTGCATCAGAATTTTTCTCT
 CACGTGGTAGACGTCACCCACCATGAAGATGCTTTAACAAAAACAAATATAATATTCGTGGCTATCCATA
 GAGAACATTACACCTCCTTGTGGGACCTGAGACATCTGCTTGTGGGCAAAATCCTCATTGATGTGAGCAA
 CAACATGAGAGTAAACCAGTACCAGAATCCAATGCAGAGTACCTGGCTTCAATTATCCCGGACTCCTTG
 ATTGTCAAAGGATTTAATGTGATCTCAGCTTGGGCACTCAGTAGGTCCTCAAGGATGCCAGCCGCCAGG
 TTTATATATGCAGCAACAATATCCAAGCTCGACAGCAGGTTATCGAGCTCGCCCGCCAGCTGAATTTTAT
 TCCTGTTGACTTGGGATCTTTGCGTCAGCCAAGGAGATTGAAAACCTTACCTCTGCGACTGTTTACTCTC
 TGAGGGGGCCAGTGGTAGTAGCCATAAGCTTGGCCACATTTTTCTTTTATTCTTTTGTGAGAGATG
 TGATACATCCATATGCCAGAAACCAGCAGAGTGACTTTTACAAGATCCCATGAGATTGTGAACAAAAC
 CTTGCCGATCGTCGCCATCACCTGCTGTCTCTGGTGTACCTGGCTGGCCTCCTGGCAGCTGCGTATCAG
 CTTTATTATGGCACTAAGTACCGCCGATTTCCCGCGTGGCTGGATACTTGGCTGCAGTGCAGGAAACAGC
 TGGGATTGCTGAGCTTCTTCTTTCAGTTGTTACAGTACGCTACAGCCTCTGCTTACCAATGAGGAGGTC
 GGAAAGATACCTGTTCTCAACATGGCTTATCAGCAGGTTATGCCAATATTGAGAACCGGTGGAACGAG
 GAGGAGGCTGGAGGATTGAGATGTACATTTCTTTGGCATCATGAGCCTGGGCTTGTGTCCCTGCTGG
 CGGTCACTTCCATCCCATCAGTGAAGCAAGCTTTGAAGTGGAGAGAGTTTCACTTCCATCCAGTCTACGCT
 TGGCTACGTCGCCCTGCTCATCAGACCTTCCAGTGTTAATTTACGGATGGAAGCGTGCCTTTGCAGAA
 GAGTACTACCGCTTTTACACACCACAAACTTCGTTCTTGCCCTCGTTTTGCCCTCCATTGTAATTTCTGG
 GTAAGATGATATTACTCTCCATGCATAAGCCGAAAGCTAAAACGAATTAAGGAGGCTGGGAAAAGAG
 CCAGTTTCTAGACGAAGGCATGGGAGGAGCGGTTCTCATCTGTCCCGAGAGAGGTCACAGTGATG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG216709 representing NM_028734
 Red=Cloning site Green=Tags(s)

MESISMMGSPKSLETFLPNGINGIKDARQVTVGVI GSGDFAKSLTIRLIRCGYHVIVIGSRNPKFASEFFP
 HVVDVTHHEDALTKTNIIFVAIHREHYTSLWDLRHLVGKILIDVSNMRVNVQYPESNAEYLAFLPDSL
 IVKGFNVI SAWALQLGPKDASRQVYICSNNIQARQVIELARQLNFI PVDLGSLSAKEIENLPLRFLTL
 WRGPVVVAISLATFFFLYSFVRDVIHPYARNQQSDFYKIPVIEIVNKTLPVAVITLLSLVYLAGLLAAAYQ
 LYYGTKYRRFPPWLDTWLQCRKQLGLLSFFFVAVHVAYSCLPMRRSERYLFLN MAYQQVHANIANEAWNE
 EEVWRIEMYISFGIMSLGLLSLLAVTSIPSVSNALN WREFSFIQSTLGYVALLITTFHVL IYGWKRAFAE
 EYRYFYTPPNFV LALVLP SIVILGKMILL LPCISRKLKRIKKGWEKSQFLDEGMGAVPHLSPERVTVM



SGP**TRRRLE** - GFP Tag - V

Restriction Sites:


Sgfl-RsrII

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_028734

ORF Size: 1467 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028734.5](#), [NP_083010.2](#)

RefSeq Size: 10396 bp

RefSeq ORF: 1470 bp

Locus ID: 74051

UniProt ID: [Q8BWB6](#)

Cytogenetics: 5 A1

Gene Summary: Metalloreductase that has the ability to reduce both Fe(3+) to Fe(2+) and Cu(2+) to Cu(1+). Uses NAD(+) as acceptor.[UniProtKB/Swiss-Prot Function]