

Product datasheet for **MG216710**

Eaf2 (NM_001113401) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eaf2 (NM_001113401) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Eaf2
Synonyms:	AW048865; Festa; FESTA-L; FESTA-S; Traits; U19
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG216710 representing NM_001113401 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGTGGACCAGCGGGACTTGCATACCTGGACCGTCGCGAGCGGGTTCTCAAGCTAGGCGAAAGTTTCG
AGAAGCAGCCGCGCTGTGCCTTCCACACCGTGCCTATGACTTCAAACCTGCTTCTATTGATACTTCTTG
TGAAGGAAATCTTGAGGTTGGCAAAGGTGAACAGGTGACAATAACTCTTCAAATATAGAAGTTCAACT
CCACCAGTCACAGTTTTCAAAGGTTCCAAGAGACCTTACTTAAAAGAATGCATTTTGATTATTAACCATG
ATACTGGGGAATGTCGCCTAGAAAAGCTCAGCAGCAACATCACTGTGAAAAACAAGAGTGAAGGGAG
TAGCAGAATCCAGTACAGACTAGAACAACAGCAACAGCAAATGTGGAATCTGCCTAGGACTTCCAATCTT
GTACAGCATTCTCCATCAGAAGAGAAGATGTCTCCAACGTCTCTAATGGATGATATTGAAAGAGAACTGA
AAGCAGAAGCTAGTCTTATGGACCAGATGAGTAGTTGTGATAGTTCATCAGATTCAAAAGTTCTTCATC
TTCAAGTAGTGAGGATAGTTCTAGTGATTCTGAAGATGATGACCAATTCTCTCTTTGGTCCAAGGAAA
TACAGCTCGGAGCACCCTAGCATGTCTGCTGGCCACAGTACAGGACTTCAGAGGCTGATGCTACTTGTC
ACCGACTTCAGGACCACAGTACCCTTCTGATGAGTACTTACGAAGTGACTTGCAGCTGAGTGAGTCAGA
AAGCGACAGTGAGGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG216710 representing NM_001113401
 Red=Cloning site Green=Tags(s)

MSGPAGLAYLDRRERVLKLGESFEKQPRCAFHTVRYDFKPASIDTSCEGNLEVGKGEQVITLPNIEGST
 PPVTVFKGSKRPYLKECILIIINHDTGECRLEKLSNITVKKTRVEGSSRIQYRLEQQQQMWNLPRTSNL
 VQHSPSEEKMSPTSLMDDIERELKAEASLMDQMSSCDSSSSDKSSSSSSSEDDDDQF SPLGPRK
 YSSEHPMSAGPQYRTSEADATCHRLQDHSTLLMSTLRSDLQLSESESDSED

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001113401

ORF Size: 786 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001113401.1 , NP_001106872.1
RefSeq Size:	1888 bp
RefSeq ORF:	789 bp
Locus ID:	106389
UniProt ID:	Q91ZD6
Cytogenetics:	16 B3
Gene Summary:	Acts as a transcriptional transactivator of ELL and ELL2 elongation activities (By similarity). Acts as a transcriptional transactivator of TCEA1 elongation activity.[UniProtKB/Swiss-Prot Function]