

## Product datasheet for **MG216755**

### Gpr19 (NM\_001167696) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr19 (NM_001167696) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gpr19
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG216755 representing NM_001167696 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGATAACGACCAGCCGCCTGTGGTACTGCCACCCTGCTGGTGCCCTTCAGAACAGCAGCTGCGCGG  
AAGCAGCTGAGGCCCTGCTGCCCATGGCTGATGGGATTGCATGAGGAGCACAGCTGGATGAGCAACAG  
GACAGAGCTTCAGTACGAGCTGAACCTGGAGAGGTGGCCACCGCCAGCATCTTCTTTGGCGCTTTGTGG  
TTGTTCTCCATCTTTGGCAATTCCTGGTATGTCTGGTCATCCACCGAGCCGGAGGACTCAGTCCACCA  
CCAATACTTCGTGGTGCATGGCGTGTGCTGACCTTTCATCAGCGTGGCCAGCAGCCGTTTGTGCGT  
GCTGCAGTTCACCACCGGAGGTGGACCTGGGTAGCGCCATGTGCAAGGTGGTCCGCTACTTCCAGTAC  
CTCACCCCGGCGTGCAGATCTACGTGCTGCTCTCCATCTGCATCGACCGCTTCTATACCATCGTCTACC  
CGCTGAGCTTCAAGGTGTCCCGAGAGAAGGCCAAGAAAATGATTGCAGCCTCCTGGATCTTGGACGCGGC  
CTTCGTGACGCCTGTCTTCTTTTCTACGGCTCTAACTGGGATAGTCACTGTAATACTTCTCCTCCCTCCC  
TCCTGGGAGGAACTGCCTATACCGTCATCCACTTCTTGGTGGGCTTCGTGATCCCTCCATCCTCATAA  
TCCTGTTTTACCAGAAGTGCATAAAGTATATCTGGAGAATAGGCACGGACGGCGGACGCTGAGAAGGAC  
GATGAACATTGTCCCAGGACAAAGGTGAAAACGGTCAAGATGTTCTGCTCTTGAACCTAGTGTTCCTG  
TTCTCCTGGCTACCTTTCATGTGGCTCAGCTCTGGCATCCCCATGAGCAAGACTACAAGAAGAGCTCCC  
TTGCTTTCACGGCAGTCAGTGGGTATCTTTCAGCTCTTCAGCCTCGAAACCCACTCTGTACTCTATTTA  
TAATGCCAATTTTCGGAGAGGGATGAAAGAGACATTCTGCATGTCTCGATGAAATGCTACCGCAGCAAT  
GCCTACACCATCAGACCAGCTCACGGATGGCCAAAAGAACTATGTGGGCATTTTCGAAAATCCCTCCC  
TGAGCAGGACCATAACCAAAGACTCCATCTATGACTCATTTGACCGGAGGCCAGGGAGAAGAAGCTCGC  
CTGGCCATCAACTCAAACCTCCAAACACTTTTGTC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


**Protein Sequence:** >MG216755 representing NM\_001167696  
 Red=Cloning site Green=Tags(s)

MDNDQPPVVTATLLVPLQNSSCAEAAEALLPHGLMGLHEEHSWMSNRTELQYELNPGEVATASIFFGALW  
 LFSIFGNLSVCLVIHRSRRTQSTTNYFVVSACADLLISVASTPFVVLQFTTGRWTLGSAMCKVRYFYQY  
 LTPGVQIYVLLSICIDRFYTIIVPLSFKVSRKAKKMAASWILDAAFVTPVFFFYGSNWDSHCNVFLPP  
 SWEGTAYTVIHFLVGFVIPSILIIILFYQKVIKIYIWRIGTDGRTLRRMTNIVPRTKVKTVMFLLNLVFL  
 FSWLPPHVAQLWHPHEQDYKSSLVFTAFTWVSFSSSASKPTLYSIYNANFRRGMKETFCMSSMKCYRSN  
 AYTIITSSRMAKRNYVGISEIPPVSRITITKDSIYDSFDREAREKKLAWPINSNPPNTFV


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001167696

<b>ORF Size:</b>	1227 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001167696.2</a>
<b>RefSeq Size:</b>	1643 bp
<b>RefSeq ORF:</b>	1230 bp
<b>Locus ID:</b>	14760
<b>UniProt ID:</b>	<a href="#">Q61121</a>
<b>Cytogenetics:</b>	6 G1
<b>Gene Summary:</b>	Orphan receptor.[UniProtKB/Swiss-Prot Function]