

Product datasheet for **MG216809**

Ldb3 (NM_001039071) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ldb3 (NM_001039071) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ldb3
Synonyms:	AW742271; PDLIM6; ZASP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG216809 representing NM_001039071
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTTACAGTGTGACTCTGACTGGGCCTGGGCCTGGGGCTCCGCTGCAGGGAGGCAAGACTTCA
 ACATGCCCTCACTATCTCCCGGATCACTCCAGGCAGCAAGGCAGCCAGTCCCAGCTCAGCCAAGGAGA
 CCTGGTGGTGGCCATTGATGGCGTCAACACAGACACCATGACCCACCTGGAGGCCAGAACAAGATCAAG
 TCGGCCAGCTACAACCTGAGCCTCACTCTACAGAAGTCCAAGCGGCCTATTCCCATCTCCACGACAGCGC
 CTCCCATCCAGTCCCGCTGCCAGTGATCCCCACCAGAAGGACCCTGCTCTGGACACGAACGGCAGCCT
 GGGGACTCCCAGCCCAGCCCTGAGGCAAGGGCCAGCCAGGGCCCTGGAGTTCGGGGACACCTTTAGC
 TCCTCCTCTCCAGACCTCTGTCTGCAGCCACTCATGGAGGCCTCGGGCCTGTCTTCTCTGGGTA
 GCCAGTGGCCAAGGCCAGCTCAGAGGGTGCACAGGGCTCAGTCAGCCCCAAAGTTCTGCCAGGCCCGAG
 CCAACCGAGGCAATATAACAACCCATTGGCCTGTACTCAGCAGAGACCCTGAGGGAGATGGCTCAGATG
 TATCAGATGAGCCTCCGAGGGAAGGCCTCAGGTGCTGGACTCCTAGGAGGGAGCCTCCCTGTGAAAGACT
 TGGCGGTGGACAGCGCCTCTCCTGTGTACCAGGCTGTGATCAAAACCCAAAGCAAGCCGGAAGATGAGGC
 TGATGAGTGGGCCCGCCGGTCTCCAACCTGCAGTCTCGCTCCTCCGCATCCTGGCCCAGATGACCGGG
 ACAGAATACATGCAAGACCCTGATGAAGAGGCTCTGCAAGGTCAAGGCCACAGGCCTCTGCTTACAGTC
 CTGCAGCCGCGCTTCTCCAGCACCTCTGCCCATACCAGCTACAGCGAGGGTCCAGTGGCCCTGCACC
 CAAGCCTCGGGTCTGAACCACTGCCAGCATCCGGCCTTCTGTCTACCAGCCAGTGCCTGCATCTTCTAC
 AGCCCATCTCCAGGAGCAATTACAGTCCAACCTCCCTATACCCCTTCGCCAGTCTCTGCTTATACCCCT
 CTCTGTCTCCACCTATACCCCTCACCTGCACCTACATACTCCCTTCCCTGCACCAGTTATACCCCT
 CTCACCTGCCCAAACCTATACCCCTACACCATCAGCTGCCTACAGTGGGGGCCCTTCAGAATCTGCCAGC
 CGTCCACCCTGGGTGACCGATGACAGCTTCTCTCAGAAGTTTGCTCCTGGAAAGGCACTACCACTGTCA
 GCAAGCAGACCCTGCCCGGGGGCTCCAGCCTACAACCCAACAGGTCCACAAGTGACACCTCTTGCCAG
 GGGCACCTTCCAGAGAGCTGAGCGCTTCCAGCTAGCAGCCGGACTCCCTCTGTGGGCACTGTAACAAT
 GTCATCAGGGGCCCTTTCTGGTGGCCATGGGCCGCTCTTGGCACCCAGAAGAGTTCAATTGTGCCTACT
 GCAAGACCTCCCTAGCAGATGTGTGCTTTGTGGAGGAGCAGAACAACGTTTACTGCGAGCGGTGTTACGA
 GCAGTTCTTTGCCCGATCTGTGCCAAATGTAACACCAAATCATGGGGAAAGTATGCACGCTCTGCGA
 CAGACATGGCACACCACCTGCTTTGTCTGTGCAGCCTGCAAGAAGCCCTTTGGGAACAGCCTCTCCACA
 TGGAGGATGGAGAGCCATCTGTGAGAAAGACTACATCAACCTGTTACGACCAAGTGCCATGGATGCGA
 TTTCCCTGTGGAAGCTGGTGACAAGTTTATCGAGGCCCTGGGGCACACCTGGCATGACACCTGCTTCATT
 TGTGCGGTCTGCCACGTGAATCTGGAGGGTCAGCCATTCTACTCTAAGAAAGACAAACCACTGTGCAAGA
 AGCACGCACATGCCATCAACGTG

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG216809 representing NM_001039071
Red=Cloning site Green=Tags(s)

```

MSYSVTLTGPWPWFRLQGGKDFNMPLTISRITPGSKAAQSQLSQGDLVVAIDGVNTDTMTHLEAQNKIK
SASYNLSLTLQKSKRPIPISTTAPPIQSPLVPIPHQKDPALDTNGSLATPSPSPEARASPGALEFGDTFS
SSFQTSVCSPLMEASGPVLPVLPVAKASSEGAQGSVSPKVLPGPSQPRQYNNP IGLYSAETLREMAQM
YQMSLRGKASGAGLLGGSLPVKDLAVDSASPVYQAVIKTQSKPEDEADEWARRSSNLQSRFRILAQMTG
TEYMQDPDEEALRRSRPQASAYSPAAAASPAASHTSYSEGPAAPKPRVTTASIRPSVYQVPVASSY
SPSPGANYSPTYTPSPAPAYTPSPAPTYTPSPAPTYTPSPAPAYTPSPAPNYTPPTSAAYSGGPSESAS
RPPWVTDDFSQKFAPGKSTTTYSKQTLPRGAPAYNPTGPQVTPPLARGTFQRAERFPASSRTPLCGHCNN
VIRGPFVAMGRSWHPEEFNCAYCKTSLADVCFVEEQNNVYCERCYEQFFAPICAKCNTKIMGEVMHALR
QTWHTTCFVCAACKKPFNGSLFHMEDGEPYCEKDYINLFSTKCHGCDFFVEAGDKFIEALGHTWHDTCFI
CAVCHVNLEGQPFYSKKDKPLCKKHAHAINV
    
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001039071

ORF Size: 1983 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039071.2](#), [NP_001034160.1](#)

RefSeq Size: 4854 bp

RefSeq ORF: 1986 bp


Locus ID: 24131

UniProt ID: [Q9JKS4](#)

Cytogenetics: 14 B

Gene Summary: May function as an adapter in striated muscle to couple protein kinase C-mediated signaling via its LIM domains to the cytoskeleton.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG216809