

Product datasheet for **MG216826**

Epsti1 (NM_178825) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Epsti1 (NM_178825) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Epsti1 |
| Synonyms: | 2310046K10Rik; 5033415K03Rik; BRES11 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >MG216826 representing NM_178825 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTACACCCGCAGTAAGGTGGTGGGCCAGGCTTGGCACTTCCTCCATCTCTCGGGACCACGCTGGAG
CTGGGCAGCGCAGGGAGCTGGGCTTGCAGCAAACCGGAGACAGAGTTTGGAGGTGGCAGCTCCCGAGGG
CCCCAAGATGGAGAGGCAGGGGCATGCTGACCAGGGAAGCGCAGGCACATACACCCTGATAGCACCAAAC
GAGAGCAGAAGGCAGAAGATACAGAGAATTGCTGAGCAAGAAGTGGCAGACCTGGAGAGATGGAAGCAAC
AGAACAGGGCAAACAGTCTACTTGGTGCCACAACGGCTGGGTGGAAGCCAGTCTGAGGCTGAAGTCAG
ACAGAAACAACAGCTGCAACAGATGCGTTCTAAATACCAGCAAAGCTCAAGAGGGATGAATCTATTCTGA
ATCAGGAAGGAAGCTGAGGAAGCAAGTTCAAAAAATGAAAGCAATCCAGAGGGAGAAGAGTAATAAAC
TAGAGGAGAAAAAGCAACTTCAAGAAGACATCAGAAGAGCCACATTTGAGAGCATCATCAGTCCAAAAAC
GGCTGAGCTCTTGAGCAGACTGGACACAGAGCGGAGAAACAGAAGTGCCTGTCTGATTGCTCCTCTGCC
ACCCAGTCTCGAGATGGGTGAGTGAAACTCACGCAGGCTGCTGTCTGCCACCTTT

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG216826 representing NM_178825
Red=Cloning site Green=Tags(s)

MYTRSKVVGPGLGTSSISRDHAGAGQRRELGLQNNRRQSLEVAAPGPKMERQGHADQGSAGTYTLIAPN
 ESRRQKIQRIAEQELADLERWKQNNRAKPVYLVPQRLGGSQSEAEVRQKQQLQQMRSKYQQKLRDESIR
 IRKEAEEAKFQMKAIQREKSNKLEEKQLQEDIRRATFREHHQSKTAELL SRLD TERRNRSACL IAPPA
 TQSSRWVSGNSRRLLSATF


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_178825

ORF Size: 687 bp

| | |
|-------------------------------|---|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_178825.4 , NP_849147.1 |
| RefSeq Size: | 1364 bp |
| RefSeq ORF: | 690 bp |
| Locus ID: | 108670 |
| UniProt ID: | Q8VDI1 |
| Cytogenetics: | 14 D3 |
| Gene Summary: | Plays a role in M1 macrophage polarization and is required for the proper regulation of gene expression during M1 versus M2 macrophage differentiation (PubMed:29217193). Might play a role in RELA/p65 and STAT1 phosphorylation and nuclear localization upon activation of macrophages (PubMed:29217193).[UniProtKB/Swiss-Prot Function] |