

Product datasheet for **MG216849**

Olfr154 (NM_013728) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr154 (NM_013728) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Olfr154
Synonyms:	912-93; MOR175-1; Olfr4-3; OR93
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG216849 representing NM_013728 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGATGCACAGAAATCAAACCTGTTGTGACTGAGTTCTTCTTCACTGGATTAACCTCCTCCTTCCATCTCC
 AAATTGTCCTCTTCTGACATTTCTCTGTGTTATCTCGCAACTCTTCTGGGGAACCTGGGAATGATCAT
 TCTTATTCACCTGGACACTCGACTCCACATCCCCATGACTTTTTCTCAGCCACTTGCCTTTGTGGAT
 GCCTGCTCCTCTTCTGTCATAAGTCCTAAGATGTTGTCTGACATGTTGTGGATAAAAAGGTGATCTCCT
 TCTTAGGTTGTGCTATCCAGCTTTGTTATTTAGCCAGTTTGTAGTGACAGAATGTTTCTCCTGGCCTC
 CATGGCTTATGATCGGTATGTTGCCATCTGTAAGCCTTTGCTGTATACTCTCATCATGTCCAGCGAGTC
 TGTGTACAGTTAGTGATTGGCCCTTACAGCATAGGCTTTATAAGCACTATGGTCCATATTATTTCTGCAT
 TTGTCCTTCCATACTGTGGTCCAAATCTCATCAATCACTTTTTCTGTGACCTTCTTCTGTCTCTCACT
 GGCATGTGCAAACTCAGATGAATAAACGTTTGCTTTTCATCGTGGCTGGGATCCTAGGTGATTCACT
 GGTATAATCATATTGGTGTCTATGTTTACATTGCCATCACCATCCTAAAGATCAGTTCTGCTGATGGGA
 GCGCAAAGCCTTCTCCACCTGCTCTTCACACCTGACAGCAGTCTCTATTCTTTACGGGACTCTCTTCTT
 TATCTATGTACGTCCAAGCTCCAGCTTCTCTGGATATCAATAAGGTTGTGCTCACTATTTTATACCACT
 GTTATCCCTATGTTGAACCCATTCAATTTATAGCCTGAGAAATAAGGAAGTCAAAGATGCATTAATCAGGA
 CATTGAAAAGCAGTTTTGCTACAGCTTGCAAGATAAAATACTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >MG216849 representing NM_013728
 Red=Cloning site Green=Tags(s)

MMHRNQTVVTEFFFTGLTSSFHLQIVLFLFCLCVYLATLLGNLGMILIHLDTRLHIPMYFFLSHLSFVD
 ACSSSVISPKMLSDFVDKKVISFLGCAIQLCLFSQFVTECFLLASMAYDRYVAICKPLLYTLIMSQRV
 CVQLVIGPYSIGFISTMVHIIISAFVLPYCGPNLINHFFCDLLPVLACANTQMNKRLLFIVAGILGVFS
 GIIILVSYVYIAITILKISSADGRRKAFSTCSSHLTAVSILYGTLFFIYVRPSSSFLDINKVVSIFYTT
 VIPMLNPFIIYSLRNKEVKDALIRTFEKQFCYSLQDKIL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_013728

ORF Size: 954 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_013728.1](#)

RefSeq Size: 1079 bp

RefSeq ORF: 957 bp


Locus ID: 27216

UniProt ID: [Q9QY00](#)

Cytogenetics: 2 D

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MG216849