

Product datasheet for **MG216852**

Kir3dl1 (NM_177749) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Kir3dl1 (NM_177749) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Kir3dl1 |
| Synonyms: | Kirl1; Kirl2; Krl1 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >MG216852 representing NM_177749 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGCTCTGGTTCCTCAGCCTCGTGTGTTCTGGTTTCTCCTAGTCCAGAGGATGTCAGCACACTTGG
GAAGCCATGATAAACCATTCTATCTGCCTGGCCAAGTTATGTTGTTCTCTGGGACAGAATGTGACTCT
TACATGTGACTCTCATAGAGGATCTAACATATTCAGCTGTACAAAGAAGAAGGGTCTCCAAACCATCAG
TTACATGAGACTACATTCAGAAGAGCCAGGTCTTTGGTCTGTGACAACAGAACATGCAGGAACCTACA
GATGCTTCCATCCTCAATATGCTAATGTATTGTCAGCACATAGTGAGCCTCTGAAAATCATAATTTCAGG
AATCTACATGAAGCCTTTCCTCTTGATTCTACAAAGTCCCTGGTGGATACAGGAGGGAATGTGACCCTG
GAATGTCACCTCAGAGAATATGTTTCGATACATACATTTTGAATTCACACAGAATGGGAATAATCAAGAACT
CTGTCCAAGTTTCTGCAGAACATCATGAAAGTGGGTACATGTCAGCTACTCAATAGGTCCCATGACACC
TGACCTTGTGGGACTTACACATGCTATGGTGCTAACAGTTACTACCCCTATGAATGGTCAGACTCCAGT
GACCCATTGACATCAAGATCACAGGTATATAAAGAAACCTTCTGTGACCCCTGATGGGCCCTGTGC
TGATGATGTCAGGAGAGACCATGACCCTGTCTTGATCTCTGACCATCAGTTTGACACGTTCCATATGTC
CAGGGAAGGGGTGCCTCGGGACAGGGGATGCCTGCAGTGCAGAGCCACAGTGGTAAATTCGAGGCAAAA
TTCCTTTAGTCCAATGATCCAAAAGGGAACATAGATGCTATGGCTTTTCAGGAACGCTTCCCATG
TGTGGTCTCTCCCAAGTACCCCTGTACCTTCTGCCAAAGGAAACTGCTGCATACACAGAAGCAGA
TACAAAACAAATAACTACAAGAACCTGCACATCTAACTGGGCTCTTAGTAACCATGGTCTTGTGCTC
ATCATCATCTTTTACTCTTGTTACTTCTCCAAACAGAATAAGTCTCAGAAACAAGCAGCAGCCAGCATGG
AACAAAGATATGAAGTAAAAACACAATAAACACACAGGACTTTGAAAGACAAGAAGGACAGGGTTTACA
TACACAGAACTTGAGCAGAGAATCTTCAACCAAACT

ACGGTACGGCGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG216852 representing NM_177749
 Red=Cloning site Green=Tags(s)

MLLWFLSLVCSGFFLVQRMSAHLGSHDKPFLSAWPSYVVPLGQNVLTLCDSHRGSNIFKLYKEEGSPNHQ
 LHETTFQKSQVFGPVTTEHAGTYRCFHPQYANVLSAHSEPLKIIISGIYMKPFLILQSPLVDTGGNVTL
 ECHSENMFDTYILISHRMGIKNSVQVSAEHHESGSHVSYSIGPMPDLVGTYTTCYGANSYYPYEWSDSS
 DPIDIKITGIYKKPSLSALMGPVLMMSGETMTLSCI SDHQFDTFHMSREGVPRGQMPAVQSHSGKFEAK
 FLLSPMIQKGNRYCYGSFRNASHVWSSPDPLYLPAKGNCPAYTEADTKTNNYKNLHILTGLLVTMVLV
 IIFYSYCYFSKQNKSKQAAASMEQEYEVKNTINTQDFERQEGQLHTQNL SRESSTKT


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_177749

| | |
|-------------------------------|--|
| ORF Size: | 1227 bp |
| OTI Disclaimer: | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p> |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_177749.4 , NP_808417.2 |
| RefSeq Size: | 1557 bp |
| RefSeq ORF: | 1230 bp |
| Locus ID: | 245616 |
| Cytogenetics: | X 58.77 cM |