

Product datasheet for **MG216854**

Sarm1 (NM_001168521) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sarm1 (NM_001168521) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Sarm1
Synonyms:	A830091115Rik; C78606; MyD885; Sarm
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG216854 representing NM_001168521
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTCTGACGCTGCTCTTCTCCGCTACAACTGTGCCGTTCTTACCATGTGAGGCCACGGCCGG
 GCGCCGATCGGCTGACAGTGCCCGACCGGATCGGAGTGGTGGCGCCAGCCCATGGTGGCTGCGGGCGG
 TCGCGGTCTCGCGAAGTGTCCGCCGAGTGGCACTGAGGTGCAAGGCCCTGGAGCGTTCGCTGCCT
 GAGCTGCAGCAGGCGCTGTCCGAGCTGAAACAGGCAAGCGCGCGGGCTGTGGCGCGGGTCTCGCCG
 AGGTCTTCCAGCTGGTAGAGGAAGCCTGGCTGCTGCCGGCGTGGCCGCGAGGTGGCCAAAGTCTATG
 CGATGCTATACGCTGGACGGTGGCCTCGACTTGTGTTGCGGCTGCTTACAGCACCGGAGCTAGAGACC
 CGTGTGCAGGCCGCGCTTGTGGAGCAGATCCTGGTGGCTGAGAACCAGGACCGCTGGCGCGCATCG
 GTCTAGGCGTGATCTTGAACCTGGCGAAGGAGCGCGAGCCTGTGAACTGGCACGAAGCGTGGCGGGCAT
 CTTGGAGCACATGTTCAAGCACTCGGAGGAGAGCTGCCAGCGGCTGGTGGCGGCCGAGGCCCTCGACCG
 GTGCTGTACTGGTCCCGCCACAGACCCGGCGCTGTGCGCCACTGCGCTCTTGCCTGGCGAAGTGGC
 CGCTGCACGGGGCCAGACGGTGCACGGTGCATGGTGGAGAAGCGCGCCGCGAGTGGCTCTTCCCGCT
 CGCTTCTCCAAGGAGGACGAGCTGCTGCGGCTGCACGCCTGCCTGGCGGTGGCGGTGTTGGCTACCAAC
 AAGGAGGTGGAACGCGAGGTGAGCATTCTGGCACATTGGCGCTTGTGAGCCGCTCGTGGCATCGCTGG
 ACCCGGCGCGTTCGCCCCTGCCTGGTGGATGCCAGTGACACAAGCCAGGGTCTGGACAGACGACCT
 GCAGAGCTGGTGTGTTGCTCGATTCTGCGCTTGGAGGCTCAGTGCATAGGAGCATTCTACCTGTGC
 GCAGAGGCTGCCATCAAGAGCCTACAGGAAAGACCAAGGTGTTACAGCGACATCGGCGCTATCCAGAGCC
 TGAACCGCTGGTTTCTACTCTACGAATGGCACACGTCGGCGCTGGCAAGCGCGCTGCGCTATT
 GGGCGAGGAGGTGCCAAGGCGCATCCTGCCCTGGTGGCCAGCTGGAAGGAAGCTGAGGTCCAGACCTGG
 CTACAGCAGATCGGCTTCTCCAGTACTGCGAGAATTTTCGGGAGCAGCAGGTAGATGGTACCTGCTTC
 TAAGACTCACAGATGAAGAACTCCAGACAGACCTAGGCATGAAATCAAGCATCACCCGCAAGAGTTCTT
 TAGGGAGCTCACAGAGCTCAAGACCTTCGCCAGCTACGCTACTTGCACCCGAGCAACCTAGCGGACTGG
 CTGGGCAGCCTGGATCCTCGTTCGCCAGTACACCTATGGCCTGGTCACTGCGGTCTGGACCGCTCCC
 TGCTGCACCGCTGTGAGAGCAGCAGCTCCTGGAGACTGTGGCATCCGCTGGGAGTGCACCGCACGG
 CATCCTCTGTGAGCCAGAGTCACTTTGCCAGACTGGCCTGAGAAGCTTGAGGAGACCAAGTCTCCAC
 GATGATGGACCCGCTGATAAGCAGTGGGAAGAGCCACCCTCACCTCCATGTCTCTTTCCTTGGCTCCAG
 AAATGCTACATTCGCCGCTGCCCTGACTGGAGGCAAGCTCAGTGGGACACCCAGATGTCTTTATCAG
 TTACCGGAGGAACTCAGGTTCCAGCTGGCCAGCCTCCTGAAGGTGCACCTGCAGCTTACAGGCTTACAG
 GTCTTTCATCGACGTGGAGAAGCTGGAAGCCGGCAAATTCGAGGACAAGCTTATCCAAAGCGTCATAGCGG
 CTCGCAATTTTGTCTGGTGTCTGTGCTGGGGCGCTGGATAAGTGCATGCAGGACCATGACTGCAAGGA
 CTGGGTGCACAAGGAGATTGTGACTGCTTAAAGCTGTGGCAAGAACATTGTGCCATCATTGATGGCTTT
 GAGTGGCTGAGCCTCAGGCGCTGCCTGAGGATATGCAGGCTGTACTCACCTTCAACGGCATCAAATGGT
 CCCATGAGTACCAGGAGGCCACCATCGAGAAGATCATCCGCTTCTACAGGGCCGCCCTCTCAGGACTC
 CTCTGCCGATCGGATACCAGTTTGGAGGGAGCTACCCAATGGGTCTGCCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG216854 representing NM_001168521
Red=Cloning site Green=Tags(s)


```
MVL TLLFSAYKLCRFF TMSGPRPGADRL TVPGPDRSGGAS PWWAAGGRGSREVSPGVGTEVQGALERSLP
ELQQALSELKQASAA RAVGAGLA E V F Q L VEEAWLL PAVGREVAQGL CDA I R L D G G L D L L L R L L Q A P E L E T
RVQAARLLEQILVAENRDRVARIGLVILNLAKEREPVELARSVAGILEHMFKHSEETCQRLVAAGGLDA
VLYWCRRTDPALLRHICALALANCALHGGQTVQRCMVEKRAAEWLFPLAFSKEDELLRLHACLAVAVLATN
KEVEREVEHSGTLALVEPLVASLDPGRFARCLVDASDTSQGRGPPDLQSLVLLLDSSRLEAQCIGAFYLC
AEAAIKSLQGKTKVFS DIGAIQSLKRLVSYSTNGTTSALAKRALRLLGEEVPRRILPCVASWKEAEVQTW
LQQIGFSQYCENFREQQVDGDL LRLTDEELQTDLGKSSITRKRFFRELTELKTFASYATCDRSNLADW
LGSLDPRFRQYTYGLVSCGLDRSLLHRVSEQQLLED CGIRLGVHRTRILSAARGHFAQTGLRSLRRPSLH
DDGPRDKQWRATLTSM SLSLAP EMLH S P L P C T G G K L S G D T P D V F I S Y R R N S G S Q L A S L L K V H L Q L H G F S
VFIDVEKLEAGKFEDKLIQSVIAARNFVLVLSAGALDKCMQDHDCKDWHVHKEIVTALSCGKNIVPIIDGF
EWPEPQALPEDMQAVLTFNGIKWSHEYQEATIEKIIIRFLQGRPSQDSSAGSDTSLEGATPMGLP
```

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja1289_f04.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:




ACCN: NM_001168521

ORF Size: 2292 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001168521.1, NP_001161993.1</u>
RefSeq Size:	5164 bp
RefSeq ORF:	2295 bp
Locus ID:	237868
UniProt ID:	<u>Q6PDS3</u>
Cytogenetics:	11 46.74 cM
Gene Summary:	Negative regulator of MYD88- and TRIF-dependent toll-like receptor signaling pathway which plays a pivotal role in activating axonal degeneration following injury. Promotes Wallerian degeneration an injury-induced axonal death pathway which involves degeneration of an axon distal to the injury site. Can activate neuronal death in response to stress. Regulates dendritic arborization through the MAPK4-JNK pathway. Involved in innate immune response. Inhibits both TICAM1/TRIF- and MYD88-dependent activation of JUN/AP-1, TRIF-dependent activation of NF-kappa-B and IRF3, and the phosphorylation of MAPK14/p38. Can restrict West Nile virus (WNV) pathogenesis.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MG216854